



## Quick-lock



## SC series



### Quick-lock cystoscopy sheath

To use with 4 mm telescope



With two fix stopcocks

Code	Length	Diameter	Notes
T6020100	224 mm	17,5 FR	With obturator ref. T6040100 and adapters ref. L4000000
T6020200	224 mm	19,5 FR	With obturator ref. T6040200 and adapters ref. L4000000
T6020300	224 mm	21 FR	With obturator ref. T6040300 and adapters ref. L4000000
T6020400	224 mm	23 FR	With obturator ref. T6040400 and adapters ref. L4000000
T6020500	224 mm	25 FR	With obturator ref. T6040500 and adapters ref. L4000000

### Quick-lock standard obturator



Code	Diameter	Notes
T6040100	17,5 FR	To use Quick-lock cystoscopy sheath ref.T6020100
T6040200	19,5 FR	To use Quick-lock cystoscopy sheath ref.T6020200
T6040300	21 FR	To use Quick-lock cystoscopy sheath ref.T6020300
T6040400	23 FR	To use Quick-lock cystoscopy sheath ref.T6020400
T6040500	25 FR	To use Quick-lock cystoscopy sheath ref.T6020500

### Quick-lock visual obturator


To use with 4 mm telescope



Code	Diameter	Notes
T6070100	17,5 FR	To use Quick-lock cystoscopy sheath ref.T6020100
T6070200	19,5 FR	To use Quick-lock cystoscopy sheath ref.T6020200
T6070300	21 FR	To use Quick-lock cystoscopy sheath ref.T6020300
T6070400	23 FR	To use Quick-lock cystoscopy sheath ref.T6020400
T6070500	25 FR	To use Quick-lock cystoscopy sheath ref.T6020500

### SC series cystoscopy sheath

To use with 4 mm telescope



With two fix stopcock

Code	Length	Diameter	Notes
T8020100	224 mm	17 FR	With obturator ref. T8040100 and adapters ref. L4000000
T8020200	224 mm	19 FR	With obturator ref. T8040200 and adapters ref. L4000000
T8020300	224 mm	20 FR	With obturator ref. T8040300 and adapters ref. L4000000
T8020400	224 mm	22 FR	With obturator ref. T8040400 and adapters ref. L4000000
T8020500	224 mm	25 FR	With obturator ref. T8040500 and adapters ref. L4000000

### SC series standard obturator



Code	Diameter	Notes
T8040100	17 FR	To use SC series cystoscopy sheath ref. T8020100
T8040200	19 FR	To use SC series cystoscopy sheath ref. T8020200
T8040300	20 FR	To use SC series cystoscopy sheath ref. T8020300
T8040400	22 FR	To use SC series cystoscopy sheath ref. T8020400
T8040500	25 FR	To use SC series cystoscopy sheath ref. T8020500

### SC series visual obturator

To use with 4 mm telescope



Code	Diameter	Notes
T8070200	19 FR	To use SC series cystoscopy sheath ref. T8070200
T8070300	20 FR	To use SC series cystoscopy sheath ref. T8070300
T8070400	22 FR	To use SC series cystoscopy sheath ref. T8070400
T8070500	25 FR	To use SC series cystoscopy sheath ref. T8070500

### Chart for instrument channel capacity - Quick-lock connecting system



Cystoscope sheaths		Capacity with albarran deflector		Capacity with adaptor bridge	
Code	Diameter	T6010010	T6010020	T6050200	T6050300
T6020100	17.5 Fr.	1 × 5 Fr.	2 × 4 Fr. or 1 × 5 Fr.	1 × 5 Fr.	2 × 4 Fr. or 1 × 5 Fr.
T6020200	19.5 Fr.	1 × 6 Fr.	2 × 5 Fr. or 1 × 6 Fr.	1 × 6 Fr.	2 × 5 Fr. or 1 × 6 Fr.
T6020300	21.0 Fr.	1 × 7 Fr.	2 × 6 Fr. or 1 × 7 Fr.	1 × 7 Fr.	2 × 6 Fr. or 1 × 7 Fr.
T6020400	23.0 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.
T6020500	25.0 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.

### Chart for instrument channel capacity - SC series connecting system




Cystoscope sheaths		Capacity with albarran deflector		Capacity with adaptor bridge	
Code	Diameter	T8010010	T8010020	T8050200	T8050300
T8020100	17.0 Fr.	1 × 5 Fr.	2 × 4 Fr. or 1 × 5 Fr.	1 × 5 Fr.	2 × 4 Fr. or 1 × 5 Fr.
T8020200	19.0 Fr.	1 × 6 Fr.	2 × 5 Fr. or 1 × 6 Fr.	1 × 6 Fr.	2 × 5 Fr. or 1 × 6 Fr.
T8020300	20.0 Fr.	1 × 7 Fr.	2 × 6 Fr. or 1 × 7 Fr.	1 × 7 Fr.	2 × 6 Fr. or 1 × 7 Fr.
T8020400	22.0 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.
T8020500	25.0 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.	1 × 9 Fr.	2 × 7 Fr. or 1 × 9 Fr.

## Quick-lock albarra deflecting mechanism

To use with 4 mm telescope




## Quick-lock adapter bridges




## SC series albarran deflecting mechanism

To use with 4 mm telescope




With one working channel  
With SC series cystoscopy sheaths  
Working channel for instrument up to 9 Fr.  
**T8010010**




With Two working channel  
With SC series cystoscopy sheaths  
Working channel for instrument up to 9 Fr.  
**T8010020**


## SC series adapter bridges



With Two working channel  
With SC series cystoscopy sheaths  
Working channel for instrument up to 9 Fr.  
**T8050300**



With one working channel  
With SC series cystoscopy sheaths  
Working channel for instrument up to 9 Fr.  
**T8050200**



Adapter bridge for diagnostic  
With SC series cystoscopy sheaths  
**T8050100**

### Quick-lock optical instrumentation

To use with 4 mm telescope



**Optical grasping forceps**  
Double action jaws for stent removal  
With Quick-lock cystoscopy sheaths form 19,5 Fr.  
**L5400035**



**Short jaw optical biopsy forceps**  
Double action jaws  
forceps located at the bottom of endoscope channel  
With Quick-lock cystoscopy sheaths form 19,5 Fr.  
**L5400007**



**Optical biopsy forceps**  
Double action jaws , spherical through cutting jaws  
With Quick-lock cystoscopy sheaths form 19,5 Fr.  
**L5400030**



**Optical scissors**  
Double action jaws  
With Quick-lock cystoscopy sheaths form 19,5 Fr.  
**L5400008**




**Optical stone crushing forceps**  
Double action jaws  
With Quick-lock cystoscopy sheaths form 23 Fr.  
**L5400006**



**Adapter**  
To use Quick-lock optical forceps  
**T5091600**

### SC series optical instrumentation

To use with 4 mm telescope



**Optical grasping forceps**  
Double action jaws for stent removal  
With SC series cystoscopy sheath form 20 Fr.  
**L8400035**



**Short jaw optical biopsy forceps**  
Double action jaws  
forceps located at the bottom of endoscope channel  
With SC series cystoscopy sheath form 20 Fr.  
**L8400031**



**Optical biopsy forceps**  
Double action jaws , spherical through cutting jaws  
With SC series cystoscopy sheath form 20 Fr.  
**L8400030**




**Optical scissors**  
Double action jaws  
With SC series cystoscopy sheath form 20 Fr.  
**L8400008**



**Optical stone crushing forceps**  
Double action jaws  
With SC series cystoscopy sheath form 25 Fr.  
**L8400006**



**Adapter**  
To use SC series optical forceps  
**T8091600**




**Adapter**  
To use SC series optical forceps  
**T8091700**

### Compatibility matrix

4 mm telescope	Optical forceps	Cystoscopy sheath	Adapter	Resectoscopy sheath
A5403000 A5403001 A5403002	L5400007 L5400008 L5400030 L5400035	T6020200 T6020300 T6020400 T6020500	—	—
		—	T5091600	T5030501 T5030521 T5040010 T5040025 T5040514 T5040518 T5040523 T5040530
	L5400006	T5020400 T5020500	—	—
		—	T5091600	T5030501 T5030521 T5040010 T5040025 T5040514 T5040518 T5040523 T5040530
	L8400008 L8400030 L8400031 L8400035	T8020300 T8020400 T8020500	—	—
		—	T8091600	T9030540 T9030550 T9040514 T9040516
		—	T8091700	T9040540 T9040550
	L8400006	T8020500	—	—
		—	T8091600	T9030550 T9040516
		—	T8091700	T9040550


### Quick-lock continuous flow laser cystoscopy



#### Contiuous flow laser cystoscopy sheath, complete

To use with 4 mm forward oblique telescope 30°, ref. A5403002  
and with 6,5Fr. working channel  
Delivered without obturator

**T6080400**



#### Quick-lock inner sheath for continuous flow laser cystoscopy

To use with 4 mm forward oblique telescope 30°, ref. A5403002  
and with 6,5Fr. working channel  
Delivered without obturator

**T6080000**



#### Quick-lock outer sheath for continuous flow laser cystoscopy

To use with quick-lock inner sheath ref. T6080000

**T6080100** 23 Fr.

Quick-lock continuous flow laser cystoscopy




Forward oblique telescope 30°, enlarged view

	Length	Diameter
A5403002	301 mm	4 mm




Quick-lock standard obturator  
 To use with inner sheath ref. T6080000  
 combi sheath ref. T6080400  
**T6080200**



Quick-lock visual obturator  
 To use with inner sheath ref. T6080000  
 combi sheath ref. T6080400  
**T6080300**

### Flexible instrumentation



#### Flexible forceps for removal of foreign body

Double action jaws, autoclavable  
length 40 cm

**L5400020** 5 Fr.  
**L5400000** 7 Fr.  
**L5400001** 9 Fr.



#### Flexible biopsy forceps

Double action jaws, autoclavable  
length 40 cm

**L5400021** 5 Fr.  
**L5400002** 7 Fr.  
**L5400003** 9 Fr.




#### Flexible scissors


Single action jaws, autoclavable  
length 40 cm

**L5400022** 5 Fr.  
**L5400004** 7 Fr.  
**L5400005** 9 Fr.

## Stone baskets




5 Fr. reusable basket for stone removal  
4 wires  
L5400100 60 cm



5 Fr. reusable basket for stone removal  
6 wires  
L5400101 60 cm

### Quick-lock overview



### SC series overview

