

Sheath and Obturator

SHEATH



OBTURATOR



Code	Length	Diameter	Notes
T1203050 T1203046	204 mm	9,0 FR	To use with 2.0 mm scope 0 deger code A3202990
T1273050 T1273046	270 mm	10,5 FR	To use with 2.7 mm scope 0 deger code A3282990
T1273052 T1273046	270 mm	10,5 FR	To use with 2.7 mm scope 30 deger code A3282992
T1293052 T1293046	270 mm	13,0 FR	To use with 2.9 mm scope 30 deger code A5293002
T1403070 T1403046	270 mm	15,0 FR	To use with 4.0 mm scope 0 deger code A3412980
T1403072 T1403046	270 mm	15,0 FR	To use with 4.0 mm scope 30 deger code A3412982

Single flow operation sheaths



With one channel and rotatable stopcock

Code	Length	Diameter	Notes
T1273150	207 mm	16,5 FR	5Fr. channel to use with 2.7 mm scope 0 deger code A3282990
T1273152	208 mm	16,5 FR	5Fr. channel to use with 2.7 mm scope 30 deger code A3282992
T1433080	207 mm	16,5 FR	5Fr. channel to use with 4.0mm scope 0 deger code A3412980
T1433102	207 mm	16,5 FR	5Fr. channel to use with 4.0mm scope 30deger code A3412982



With two channel and one rotatable stopcock

Code	Length	Diameter	Notes
T1403030	215 mm	22,5 FR	7Fr. channel to use with 4.0 mm scope 0 deger code A3412980
T1403052	215 mm	22,5 FR	7Fr. channel to use with 4.0 mm scope 30 deger code A3412982

Continuous flow operation sheaths



With two fixed stopcock

Code	Length	Diameter	Notes
T1202470	148 mm	15,0 FR	5Fr. channel to use with 2.0 mm scope 0 deger code A3202990
T1403100	217 mm	22,5 FR	7Fr. channel to use with 4.0 mm scope 0 deger code A3412980
T1403102	215 mm	22,5 FR	7Fr. channel to use with 4.0mm scope 30 deger code A3412982



With two rotatable stopcock

Code	Length	Diameter	Notes
T1273202	220 mm	14,9 FR	5Fr. channel to use with 2.7 mm scope 30 deger code A3282992

Albarran deflecting mechanism



Quick-lock hysteroscopy outer sheath
single flow with one rotatable stopcock

Code	Length	Diameter	Notes
T602030I	182 mm	21 FR	To use with Albarran deflecting mechanism ref.T601001I



Standard obturator quick-lock


Code	Diameter	Notes
T604030I	21 FR	To use with hysteroscopy outer sheath ref.T602030I



Albarran deflecting mechanism quick-lock

Code	Diameter	Notes
T601001I	7 FR	To use with 4 mm telescope and hysteroscopy outer sheath ref.T602030I

Instrumentation



Flexible forceps for removal of foreign body

Double action jaws, autoclavable
length 40 cm

L5400020 5 Fr.

L5400000 7 Fr.



Flexible biopsy forceps

Double action jaws, autoclavable
length 40 cm

L5400021 5 Fr.

L5400002 7 Fr.



Flexible scissors

Single action jaws, autoclavable
length 40 cm

L5400022 5 Fr.

L5400004 7 Fr.



Flexible IUD grasping forceps

Double action jaws
length 40 cm

L0700400 5 Fr.

L0700410 7 Fr.

HF electrodes



Monopolar button electrode

length 45cm

L5400046 3 Fr.

L5400047 4 Fr.

L5400009 5 Fr.

L5400010 6 Fr.

L5400011 7 Fr.



Monopolar needle electrode

length 45cm

L5400048 3 Fr.

L5400013 5 Fr.

L5400014 7 Fr.



Monopolar loop electrode


length 45cm

L5400017 5 Fr.

L5400015 7 Fr.

Note : Identify your electrosurgical generator brand and required cable length 3 meter or 5 meter .

Laparoscopic ring applicators



For single puncture approach

Code	Length	Diameter	Notes
T1060331	33 mm	6 mm	To use with operative laparoscope with 6 mm working channel




For double puncture approach

Code	Length	Diameter	Notes
T1060341	30 mm	6,5 mm	To use with 6,5 mm trocars sleeve ref. T1071070


Adapters



SC series adapter
for record syringe
T9091000



Quick-lock adapter
with safety lock mechanism
T5091001



Luer-lock connector
with stopcock, dismantling
L4000001



Luer-lock connector, male
Tube diameter 9mm
L4000000




Quick-Lock adapter
for record syringe
T9091000




Catheter adapter without
stopcock, Luer-lock
L4000002

Bougies



Bougies a boule
9-12-15-18 Fr.
L5410004



Bougies a boule
21-24-27-30 Fr.
L5410005

Uterine manipulator



“COHEN” intra uterin ptobe
With two adaptors
L1060700

Uterine manipulator



Uterine manipulator handle

L5005090

Complete uterine manipulator
including six working inserts.

L5005100



Working insert long with conus

L5005092



Working insert short with conus

L5005091



Working insert 3×45 mm

L5005094



Working insert 3×55 mm

L5005093



Working insert 8×80 mm

L5005096



Working insert 10×100 mm

L5005095