


General Information





Code	Diameter	Length	Angel	Note	* Recmadad use
A3202600	2,0 mm	260 mm	0	WA	G, D, H, L
A3202990	2,0 mm	235 mm	0		G, D, H, L
A3282990	2,7 mm	301 mm	0	WA	G, D, H
A3282992	2,7 mm	301 mm	30		G, D, H
A5293000	2,9 mm	301 mm	0		G, U
A5293001	2,9 mm	301 mm	12		G, U
A5293002	2,9 mm	301 mm	30		G, U, D, H, L
A5303001	3,0 mm	301 mm	12		G, U
A5303002	3,0 mm	301 mm	30		G, U
A3312990	3,0 mm	301 mm	0	WA	G, L
A3312992	3,0 mm	301 mm	30	WA	G, L
A3412980	4,0 mm	301 mm	0		G, D, H, L
A3412981	4,0 mm	301 mm	12		G
A3412982	4,0 mm	301 mm	30		G, D, H
A3412990	4,0 mm	301 mm	0	WA	G
A3412992	4,0 mm	301 mm	30	WA	G
A2053040	5,0 mm	306 mm	0		G, L
A2053042	5,0 mm	306 mm	30		G, L
A2053043	5,0 mm	306 mm	45		L
A2053150	5,5 mm	316 mm	0		G, L
A2053152	5,5 mm	316 mm	30		G, L
A2053153	5,5 mm	316 mm	45		G, L
A2054500	5,5 mm	458 mm	0		G, L
A2054502	5,5 mm	458 mm	30		G, L
A2054503	5,5 mm	458 mm	45		G, L
A2103250	10,0 mm	325 mm	0	HD	G, L
A2103252	10,0 mm	325 mm	30	HD	G, L
A2103253	10,0 mm	325 mm	45	HD	G, L
A2104500	10,0 mm	450 mm	0		G, L
A2104502	10,0 mm	450 mm	30		G, L
A2104503	10,0 mm	450 mm	45		G, L
A5211870	1,9 / 2,1 mm	186 mm	0		U
A5271870	2,7 mm	186 mm	0		U
A5271872	2,7 mm	187 mm	30		U
A5271875	2,7 mm	188 mm	70		U
A5403000	4,0 mm	301 mm	0		U
A5403002	4,0 mm	301 mm	30		U
A5403001	4,0 mm	301 mm	12		U
A5403005	4,0 mm	303 mm	70		U

* WA = Wide Angel

* HD = High definition

* G = Gynecology

* L = Laproscopy

* U = Urology

* H = Operative hysteroscopy


* D = Diagnostc



Code	Diameter	Length	Angel	With working channel	* Recmaded use
A2102730	10,0 mm	273 mm	0	5,0	G, L
A2102740	10,0 mm	273 mm	0	6,0	G, L
A2122771	12,0 mm	273 mm	12	7,5	G, L


* L = Laproscopy

* G = Gynecology



Code	Diameter	Length	Angel	With working channel	* Recmaded use
A2102760	10,0 mm	273 mm	0	6,0	G, L

* L = Laproscopy
* G = Gynecology



Code	Diameter	Length	Angel	With working channel	* Recmaded use
A5800010	8,5mm × 6,0 mm	222 mm	0	13 Fr	H
A5800000	8,3mm × 6,0 mm	222 mm	0	13 Fr	U
A5800002	8,3mm × 6,0 mm	226 mm	25	13 Fr	U
A5802630	8,3mm × 5,25 mm	264 mm	0	15 Fr	U

* U = Urology

* H = Operative hysteroscopy