


SPECIALITY LARYNGOSCOPE BLADES

MACINTOSH POLIO


A modification of the English Macintosh blade which is offset from the handle at an obtuse angle to allow intubation of patients on respirators, in body jackets, and other difficult situations. It is a curved blade type for indirect laryngoscope.



Blade Ref.	Size	A	B
H-251-03	3	130 mm	108 mm

OXFORD

Suitable for premature infants, babies, and children up to the age of four. There is sufficient overhang on the open side to prevent the lips from obscuring vision and the broad, flat lower surface is a help in the small child with an extreme degree of cleft palate.




Blade Ref.	Size	A	B
H-261-01	I	102 mm	80 mm

SPECIALITY LARYNGOSCOPE BLADES

MILLER PORT

The Miller-Port is a modification of the Miller blade that allows the individual to perform laryngoscopy and intubation with less risk of hypoxia. Built into the Miller 0 and Miller I blade is a tube that permits delivery of oxygen or other gas mixtures during intubation.


- Autoclavable up to 134°C approximately 4,000 times. No adverse affect in using ETO, Steam Autoclaving.



Blade Ref.	Size	A	B
H-291-00	0 (Neonate)	78 mm	53 mm
H-291-01	I (Infant)	102 mm	80 mm

SUN-FLEX FLEXIBLE TIP LARYNGOSCOPE BLADES


- Available in size 2,3 & 4.
- Movable tip facilitates difficult intubation.
- Reduced likelihood of damage to teeth and soft tissues.
- Autoclavable up to 134°C approximately 4,000 times. No adverse affect in using ETO, Steam Autoclaving.
- An enhanced viewing area compared to other brands due to channel in the lower flange.



Blade Ref.	Size	A	B
H-301-02	2	100 mm	78 mm
H-301-03	3	130 mm	108 mm
H-301-04	4	155 mm	133 mm

STANDARD 'D' LED MACINTOSH BLADES



- Suitable for everyday hospital use
- LED provides bright illumination upto 4000 Lux
- Available in size 1, 2, 3 & 4.
- Made from Austenitic Stainless steel thereby making them rust proof.
- Perfect solution to cleaning problems confronting EMS field



Blade Ref.	Size	Description	A	B
H-441-01	1	Standard 'D' Mac Blade	90 mm	66 mm
H-441-02	2	Standard 'D' Mac Blade	100 mm	76 mm
H-441-03	3	Standard 'D' Mac Blade	130 mm	106 mm
H-441-04	4	Standard 'D' Mac Blade	153 mm	129 mm

STANDARD 'D' LED MILLER BLADES

- Suitable for everyday hospital use.
- LED provides bright illumination upto 4000 Lux.
- Available in size 00, 0,1,2,3 & 4.
- Made from Austenitic Stainless steel thereby making them rust proof.
- Perfect solution to cleaning problems confronting EMS field.



Blade Ref.	Size	Description	A	B
H-461-00	0	Standard 'D' Miller Blade	79mm	54mm
H-461-01	1	Standard 'D' Miller Blade	105mm	80mm
H-461-02	2	Standard 'D' Miller Blade	157mm	133mm
H-461-03	3	Standard 'D' Miller Blade	197mm	172mm
H-461-04	4	Standard 'D' Miller Blade	208mm	184mm
H-461-05	00	Standard 'D' Miller Blade	65mm	41mm

MACINTOSH F.O. LARYNGOSCOPE SETS

Ref.	Set Detail	Case Ref.
H-801-01	Macintosh Set of 3, Consisting of Blade 1,2,3 with 2.5V Xenon Battery Handle, Medium	C-084-01
H-801-02	Macintosh Set of 3, Consisting of Blade 2,3,4 with 2.5V Xenon Battery Handle, Medium	C-084-01
H-801-03	Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 2.5V Xenon Battery Handle, Medium	C-084-01
H-891-01	Macintosh Set of 3, Consisting of Blade 1,2,3 with 2.5V LED Battery Handle, Medium	C-084-01
H-891-02	Macintosh Set of 3, Consisting of Blade 2,3,4 with 2.5V LED Battery Handle, Medium	C-084-01
H-891-03	Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 2.5V LED Battery Handle, Medium	C-084-01
H-801-11	Macintosh Set of 3, Consisting of Blade 1,2,3 with 3.5V Xenon Battery Handle, Medium	C-084-01
H-801-12	Macintosh Set of 3, Consisting of Blade 2,3,4 with 3.5V Xenon Battery Handle, Medium	C-084-01
H-801-13	Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 3.5V Xenon Battery Handle, Medium	C-084-01
H-891-11	Macintosh Set of 3, Consisting of Blade 1,2,3 with 3.5V LED Battery Handle, Medium	C-084-01
H-891-12	Macintosh Set of 3, Consisting of Blade 2,3,4 with 3.5V LED Battery Handle, Medium	C-084-01
H-891-13	Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 3.5V LED Battery Handle, Medium	C-084-01



MILLER F.O. LARYNGOSCOPE SETS

Ref.	Set Detail	Case Ref.
H-811-01	Miller Set of 3, Consisting of Blade 1,2,3 with 2.5V Xenon Battery Handle, Medium	C-084-02
H-811-02	Miller Set of 3, Consisting of Blade 00,0,1 with 2.5V Xenon Battery Handle, Small	C-084-03
H-811-03	Miller Set of 4, Consisting of Blade 1,2,3,4 with 2.5V Xenon Battery Handle, Medium	C-084-02
H-811-04	Pediatric Set of 3 of Blades Miller 00,0 & Macintosh 0 with 2.5V Xenon Battery Handle, Small	C-084-03
H-901-01	Miller Set of 3, Consisting of Blade 1,2,3 with 2.5V LED Battery Handle, Medium	C-084-02
H-901-02	Miller Set of 3, Consisting of Blade 00,0,1 with 2.5V LED Battery Handle, Small	C-084-03
H-901-03	Miller Set of 4, Consisting of Blade 1,2,3,4 with 2.5V LED Battery Handle, Medium	C-084-02
H-901-04	Pediatric Set of 3 of Blades Miller 00,0 & Macintosh 0 with 2.5V LED Battery Handle, Small	C-084-03
H-811-11	Miller Set of 3, Consisting of Blade 1,2,3 with 3.5V Xenon Battery Handle, Medium	C-084-02
H-811-12	Miller Set of 3, Consisting of Blade 00,0,1 with 3.5V Xenon Battery Handle, Small	C-084-03
H-811-13	Miller Set of 4, Consisting of Blade 1,2,3,4 with 3.5V Xenon Battery Handle, Medium	C-084-02
H-811-14	Pediatric Set of 3 of Blades Miller 00,0 & Macintosh 0 with 3.5V Xenon Battery Handle, Small	C-084-03
H-901-11	Miller Set of 3, Consisting of Blade 1,2,3 with 3.5V LED Battery Handle, Medium	C-084-02
H-901-12	Miller Set of 3, Consisting of Blade 00,0,1 with 3.5V LED Battery Handle, Small	C-084-03
H-901-13	Miller Set of 4, Consisting of Blade 1,2,3,4 with 3.5V LED Battery Handle, Medium	C-084-02
H-901-14	Pediatric Set of 3 of Blades Miller 00,0 & Macintosh 0 with 3.5V LED Battery Handle, Small	C-084-03



GREENLINE SUNFLEX F.O. LARYNGOSCOPE SETS

Ref.	Set Detail	Case Ref.
H-821-01	Sunflex Set of 1 Blade, Macintosh 3 with 2.5V Xenon Battery Handle, Medium	C-084-04
H-821-02	Sunflex Set of 1 Blade, Macintosh 4 with 2.5V Xenon Battery Handle, Medium	C-084-04
H-821-03	Sunflex Set of 2 Blades, Macintosh 3,4 with 2.5V Xenon Battery Handle, Medium	C-084-04
H-911-01	Sunflex Set of 1 Blade, Macintosh 3 with 2.5V LED Battery Handle, Medium	C-084-04
H-911-02	Sunflex Set of 1 Blade, Macintosh 4 with 2.5V LED Battery Handle, Medium	C-084-04
H-911-03	Sunflex Set of 2 Blades, Macintosh 3,4 with 2.5V LED Battery Handle, Medium	C-084-04
H-821-11	Sunflex Set of 1 Blade, Macintosh 3 with 3.5V Xenon Battery Handle, Medium	C-084-04
H-821-12	Sunflex Set of 1 Blade, Macintosh 4 with 3.5V Xenon Battery Handle, Medium	C-084-04
H-821-13	Sunflex Set of 2 Blades, Macintosh 3,4 with 3.5V Xenon Battery Handle, Medium	C-084-04
H-911-11	Sunflex Set of 1 Blade, Macintosh 3 with 3.5V LED Battery Handle, Medium	C-084-04
H-911-12	Sunflex Set of 1 Blade, Macintosh 4 with 3.5V LED Battery Handle, Medium	C-084-04
H-911-13	Sunflex Set of 2 Blades, Macintosh 3,4 with 3.5V LED Battery Handle, Medium	C-084-04



MACINTOSH LARYNGOSCOPE SETS WITH INTEGRATED FIBERS

Ref.	Set Detail	Case Ref.
H-861-01	Maxlite Macintosh Set of 3, Consisting of Blade 1,2,3 with 2.5V Xenon Battery Handle, Medium	C-088-01
H-861-02	Maxlite Macintosh Set of 3, Consisting of Blade 2,3,4 with 2.5V Xenon Battery Handle, Medium	C-088-01
H-861-03	Maxlite Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 2.5V Xenon Battery Handle, Medium	C-088-01
H-921-01	Maxlite Macintosh Set of 3, Consisting of Blade 1,2,3 with 2.5V LED Battery Handle, Medium	C-088-01
H-921-02	Maxlite Macintosh Set of 3, Consisting of Blade 2,3,4 with 2.5V LED Battery Handle, Medium	C-088-01
H-921-03	Maxlite Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 2.5V LED Battery Handle, Medium	C-088-01
H-861-11	Maxlite Macintosh Set of 3, Consisting of Blade 1,2,3 with 3.5V Xenon Battery Handle, Medium	C-088-01
H-861-12	Maxlite Macintosh Set of 3, Consisting of Blade 2,3,4 with 3.5V Xenon Battery Handle, Medium	C-088-01
H-861-13	Maxlite Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 3.5V Xenon Battery Handle, Medium	C-088-01
H-921-11	Maxlite Macintosh Set of 3, Consisting of Blade 1,2,3 with 3.5V LED Battery Handle, Medium	C-088-01
H-921-12	Maxlite Macintosh Set of 3, Consisting of Blade 2,3,4 with 3.5V LED Battery Handle, Medium	C-088-01
H-921-13	Maxlite Macintosh Set of 4, Consisting of Blade 1,2,3,4 with 3.5V LED Battery Handle, Medium	C-088-01



MILLER LARYNGOSCOPE SETS WITH INTEGRATED FIBERS

Ref.	Set Detail	Case Ref.
H-871-01	Maxlite Miller Set of 3, Consisting of Blade 1,2,3 with 2.5V Xenon Battery Handle, Medium	C-088-02
H-871-02	Maxlite Miller Set of 3, Consisting of Blade 00,0,1 with 2.5V Xenon Battery Handle, Small	C-088-03
H-871-03	Maxlite Miller Set of 4, Consisting of Blade 1,2,3,4 with 2.5V Xenon Battery Handle, Medium	C-088-02
H-931-01	Maxlite Miller Set of 3, Consisting of Blade 1,2,3 with 2.5V LED Battery Handle, Medium	C-088-02
H-931-02	Maxlite Miller Set of 3, Consisting of Blade 00,0,1 with 2.5V LED Battery Handle, Small	C-088-03
H-931-03	Maxlite Miller Set of 4, Consisting of Blade 1,2,3,4 with 2.5V LED Battery Handle, Medium	C-088-02
H-871-11	Maxlite Miller Set of 3, Consisting of Blade 1,2,3 with 3.5V Xenon Battery Handle, Medium	C-088-02
H-871-12	Maxlite Miller Set of 3, Consisting of Blade 00,0,1 with 3.5V Xenon Battery Handle, Small	C-088-03
H-871-13	Maxlite Miller Set of 4, Consisting of Blade 1,2,3,4 with 3.5V Xenon Battery Handle, Medium	C-088-02
H-931-11	Maxlite Miller Set of 3, Consisting of Blade 1,2,3 with 3.5V LED Battery Handle, Medium	C-088-02
H-931-12	Maxlite Miller Set of 3, Consisting of Blade 00,0,1 with 3.5V LED Battery Handle, Small	C-088-03
H-931-13	Maxlite Miller Set of 4, Consisting of Blade 1,2,3,4 with 3.5V LED Battery Handle, Medium	C-088-02



FLEXITIP LARYNGOSCOPE SETS WITH INTEGRATED FIBERS

Ref.	Set Detail	Case Ref.
H-881-01	Maxlite Flexitip Set of 1 Blade, Macintosh 3 with 2.5V Xenon Battery Handle, Medium	C-088-04
H-881-02	Maxlite Flexitip Set of 1 Blade, Macintosh 4 with 2.5V Xenon Battery Handle, Medium	C-088-04
H-881-03	Maxlite Flexitip Set of 2 Blades, Macintosh 3,4 with 2.5V Xenon Battery Handle, Medium	C-088-04
H-941-01	Maxlite Flexitip Set of 1 Blade, Macintosh 3 with 2.5V LED Battery Handle, Medium	C-088-04
H-941-02	Maxlite Flexitip Set of 1 Blade, Macintosh 4 with 2.5V LED Battery Handle, Medium	C-088-04
H-941-03	Maxlite Flexitip Set of 2 Blades, Macintosh 3,4 with 2.5V LED Battery Handle, Medium	C-088-04
H-881-11	Maxlite Flexitip Set of 1 Blade, Macintosh 3 with 3.5V Xenon Battery Handle, Medium	C-088-04
H-881-12	Maxlite Flexitip Set of 1 Blade, Macintosh 4 with 3.5V Xenon Battery Handle, Medium	C-088-04
H-881-13	Maxlite Flexitip Set of 2 Blades, Macintosh 3,4 with 3.5V Xenon Battery Handle, Medium	C-088-04
H-941-11	Maxlite Flexitip Set of 1 Blade, Macintosh 3 with 3.5V LED Battery Handle, Medium	C-088-04
H-941-12	Maxlite Flexitip Set of 1 Blade, Macintosh 4 with 3.5V LED Battery Handle, Medium	C-088-04
H-941-13	Maxlite Flexitip Set of 2 Blades, Macintosh 3,4 with 3.5V LED Battery Handle, Medium	C-088-04



MACINTOSH STANDARD LARYNGOSCOPE SETS



Ref.	Set Detail	Case Ref.
H-831-01	Macintosh Set of 3, Consisting of Blade 1, 2, 3 and Battery Handle Medium	C-086-01
H-831-02	Macintosh Set of 3, Consisting of Blade 2, 3, 4 and Battery Handle Medium	C-086-01
H-831-03	Macintosh Set of 4, Consisting of Blade 1, 2, 3, 4 and Battery Handle Medium	C-086-01

MILLER STANDARD LARYNGOSCOPE SETS



Ref.	Set Detail	Case Ref.
H-841-01	Miller Set of 3 Consisting of Blade 1, 2, 3 and Battery Handle Medium	C-086-02
H-841-02	Miller Set of 3 Consisting of Blade 00,0,1 and Battery Handle Small	C-086-03
H-841-03	Miller Set of 4 Consisting of Blade 1, 2, 3, 4 and Battery Handle Medium	C-086-02
H-841-04	Pediatric Set of 3 Blades, Miller 00, 0 Macintosh 0 and Battery Handle small	C-086-03

SUNFLEX STANDARD LARYNGOSCOPE SETS



Ref.	Set Detail	Case Ref.
H-851-01	Sunflex Set of 1 Blade, Macintosh Blade 3 and Battery Handle Medium	C-086-04
H-851-02	Sunflex Set of 1 Blade, Macintosh Blade 4 and Battery Handle Medium	C-086-04
H-851-03	Sunflex Set of 2 Blades, Macintosh Blade 3, 4 and Battery Handle Medium	C-086-04

STANDARD DIAGNOSTIC SETS



Opaline / Econoline Otoscopes and Optek Ophthalmoscopes are equipped with an easy to use Bayonet locking / Screw Locking System, adaptable to same Handle for 2.5V and 3.5V vacuum light, fixed with a rheostat to adjust the intensity of light.

The power sources are designed to be sturdy and can be powered by regular batteries as well as 2.5V / 3.5V Rechargeable Batteries without needing additional components.

The Knurled Handles made surface a secure grip in the Hand, simple changes of batteries at the Base of the Handles.

STANDARD F.O. OTOSCOPE SET

Complete with Standard F.O. Otoscope, Standard Battery Handle, 1 Set of 5 Reusable Eartip 2.5, 3.5, 4.5, 5.5 & 9.0mm, 7 each 2.5mm & 4.0mm Disposable Eartips with Spare Lamp, packed in Case.

Description	Ref (2.5V)	Ref (3.5V)
Standard F.O. Otoscope Set, 2.5V	O-600-25	
Standard F.O. Otoscope Set, 3.5V, with 3.5V Li-ion Rechargeable Battery and Zoom Charging Station		O-601-35
Standard F.O. Otoscope Set, 3.5V, with Mains plug in handle including 3.5V Li-ion Rechargeable Battery		O-602-35
Standard F.O. Otoscope Set, 3.5V, W/O Zoom Charging Station		O-603-35
Replacement Lamp		
Bulb for Otoscope 2.5V, Xenon/Halogen	S-033-25	
Bulb for Otoscope 3.5V, Xenon/Halogen		S-070-35



STANDARD F.O. OPHTHALMOSCOPE SET

Complete with Standard F.O. Ophthalmoscope, Standard Battery Handle with Spare Lamp, packed in Case.



Description	Ref (2.5V)	Ref (3.5V)
Standard F.O. Ophthalmoscope Set, 2.5V	O-605-25	
Standard F.O. Ophthalmoscope Set, 3.5V, with 3.5V Li-ion Rechargeable Battery Medium 3.5V, and Zoom Charging Station		O-606-35
Standard F.O. Ophthalmoscope Set, 3.5V, with Mains plug in handle including 3.5V Li-ion Rechargeable Battery		O-607-35
Standard F.O. Ophthalmoscope Set, 3.5V, W/O Zoom Charging Station		O-608-35
Replacement Lamp		
Bulb for Ophthalmoscope 2.5V	S-029-25	
Bulb for Ophthalmoscope 3.5V		S-031-35

* Batteries are not included in the Sets.

STANDARD F.O. OTO / OPHTHALMOSCOPE SET

Complete with Standard F.O. Oscope, Standard F.O. Ophthalmoscope, Standard Battery Handle, 1 Set of 5 Reusable Eartip 2.5, 3.5, 4.5, 5.5 & 9.0mm, 7 each 2.5mm & 4.0mm Disposable Eartips with Spare Lamp, packed in Case.

Description	Ref (2.5V)	Ref (3.5V)
Standard F.O. Oto / Ophthalmoscope Set, 2.5V	O-610-25	
Standard F.O. Oto / Ophthalmoscope Set, 3.5V with 3.5V Li-ion Rechargeable Battery and Zoom Charging Station		O-611-35
Standard F.O. Oto / Ophthalmoscope Set, 3.5V, with Mains plug-in handle including 3.5V Li-ion Rechargeable Battery		O-612-35
Standard F.O. Oto/Ophthalmoscope Set, 3.5V, W/O Zoom Charging Station		O-613-35
Replacement Lamp		
Bulb for Oscope 2.5V, Xenon/Halogen	S-033-25	
Bulb for Oscope 3.5V, Xenon/Halogen		S-070-35
Bulb for Ophthalmoscope 2.5V	S-029-25	
Bulb for Ophthalmoscope 3.5V		S-031-35



SPECIAL F.O. OTO / OPHTHALMOSCOPE SET

Complete with Special F.O. Oscope, Special F.O. Ophthalmoscope. Two Special Battery Handles, Special Nasal Illuminator with F.O. Light, 1 Set of 5 Reusable Eartip 2.5, 3.5, 4.5, 5.5 & 9.0mm, 7 each 2.5mm & 4.0mm Disposable Eartips with Spare Lamp, packed in Case.



Description	Ref (2.5V)	Ref (3.5V)
Special F.O. Oto / Ophthalmoscope Set, 2.5V	O-615-25	
Special F.O. Oto / Ophthalmoscope Set, 3.5V, with 3.5V Li-ion Rechargeable Battery and Zoom Charging Station		O-616-35
Special F.O. Oto / Ophthalmoscope Set, 3.5V, with Mains plug-in handle including 3.5V Li-ion Rechargeable Battery		O-617-35
Special F.O. Oto/Ophthalmoscope Set, 3.5V, W/O Zoom Charging Station		O-618-35
Replacement Lamp		
Bulb for Oscope 2.5V, Xenon/Halogen	S-033-25	
Bulb for Oscope 3.5V, Xenon/Halogen		S-070-35
Bulb for Ophthalmoscope 2.5V	S-029-25	
Bulb for Ophthalmoscope 3.5V		S-031-35

* Batteries are not included in the Sets.

ECONOMY OTO / OPHTHALMOSCOPE SET

Complete with Otoscope, Ophthalmoscope, 1 Pocket Handle, 7 each Disposable Specula Dia 2.5mm & 4.0mm, in Plastic Case.



Description	Ref
Complete Set with Spare Lamps	O-460-01
Complete Set without Spare Lamps	O-460-02
Replacement Lamp	
For Otoscope	S-033-25
For Ophthalmoscope	S-059-25

DELUXE OTO / OPHTHALMOSCOPE SET

Complete with Otoscope, Ophthalmoscope, 2 Pocket Handles, 7 each Disposable Specula Dia 2.5mm & 4.0mm, in Plastic Case / Zipper Pouch.



Description	Ref
Complete Set with Spare Lamps	O-458-01
Complete Set without Spare Lamps	O-458-02
Complete Set with Spare Lamp in Pouch	O-462-01
Complete Set without Spare Lamp in Pouch	O-462-02
Replacement Lamp	
For Otoscope	S-033-25
For Ophthalmoscope	S-059-25

* Batteries are not included in the Sets.