



NX SYSTEM

NX PRIMARY KNEE (PS & CR)

NX Primary Knee(PS&CR)





NX-PS Femoral Condyle

Code	Description	Material
5210-01011	PS Femoral Condyle, 1L#,width 57.5mm	CoCrMo
5210-01012	PS Femoral Condyle, 1R#,width 57.5mm	
5210-01021	PS Femoral Condyle, 2L#, width 58.5mm	
5210-01022	PS Femoral Condyle, 2R#, width 58.5mm	
5210-01031	PS Femoral Condyle, 3L#, width 61mm	
5210-01032	PS Femoral Condyle, 3R#, width 61mm	
5210-01041	PS Femoral Condyle, 4L#, width 63mm	
5210-01042	PS Femoral Condyle, 4R#, width 63mm	
5210-01051	PS Femoral Condyle, 5L#, width 67.5mm	
5210-01052	PS Femoral Condyle, 5R#, width 67.5mm	
5210-01061	PS Femoral Condyle, 6L#, width 68.5mm	
5210-01062	PS Femoral Condyle, 6R#, width 68.5mm	
5210-01071	PS Femoral Condyle, 7L#, width 71mm	
5210-01072	PS Femoral Condyle, 7R#, width 71mm	
5210-01081	PS Femoral Condyle, 8L#, width 73mm	
5210-01082	PS Femoral Condyle, 8R#, width 73mm	



NX-CR Femoral Condyle

Code	Description	Material
5214-01011	CR Femoral Condyle, 1L#,width 57.5mm	CoCrMo
5214-01012	CR Femoral Condyle,1R#,width 57.5mm	
5214-01021	CR Femoral Condyle, 2L#,width 58.5mm	
5214-01022	CR Femoral Condyle,2R#, width 58.5mm	
5214-01031	CR Femoral Condyle,3L#, width 61mm	
5214-01032	CR Femoral Condyle, 3R#, width 61mm	
5214-01041	CR Femoral Condyle, 4L#, width 63mm	
5214-01042	CR Femoral Condyle, 4R#, width 63mm	
5214-01051	CR Femoral Condyle, 5L#, width 67.5mm	
5214-01052	CR Femoral Condyle, 5R#, width 67.5mm	
5214-01061	CR Femoral Condyle, 6L#, width 68.5mm	
5214-01062	CR Femoral Condyle, 6R#, width 68.5mm	
5214-01071	CR Femoral Condyle, 7L#, width 71mm	
5214-01072	CR Femoral Condyle, 7R#, width 71mm	
5214-01081	CR Femoral Condyle, 8L#, width 73mm	
5214-01082	CR Femoral Condyle, 8R#, width 73mm	



NX-PS Insert

Code	Description	Material
5217-10011	PS Tibial Insert, 1#/2#, width 61mm ,height 10mm	UHMWPE
5217-12011	PS Tibial insert, 1#/2#, width 61mm,height 12mm	
5217-14011	PS Tibial insert, 1#/2#, width 61mm,height 14mm	
5217-16011	PS Tibial insert, 1#/2#, width 61mm,height 16mm	
5217-18011	PS Tibial insert, 1#/2#, width 61mm,height 18mm	
5217-10021	PS Tibial insert, 3#/4#,width 66mm, height 10mm	
5217-12021	PS Tibial insert, 3#/4#, width 66mm,height 12mm	
5217-14021	PS Tibial insert, 3#/4#,width 66mm, height 14mm	
5217-16021	PS Tibial insert, 3#/4#, width 66mm,height 16mm	
5217-18021	PS Tibial insert,3#/4#,width 66mm, height 18mm	
5217-10031	PS Tibial insert, 5#/6#, width 70mm,height 10mm	
5217-12031	PS Tibial insert,5#/6#, width 70mm,height 12mm	
5217-14031	PS Tibial insert, 5#/6#, width 70mm,height 14mm	
5217-16031	PS Tibial insert, 5#/6#, width 70mm,height 16mm	
5217-18031	PS Tibial insert,5#/6#, width 70mm,height 18mm	
5217-10041	PS Tibial insert, 7#/8#, width 75mm,height 10mm	
5217-12041	PS Tibial insert, 7#/8#, width 75mm,height 12mm	
5217-14041	PS Tibial insert, 7#/8#, width 75mm,height 14mm	
5217-16041	PS Tibial insert, 7#/8#, width 75mm,height 16mm	
5217-18041	PS Tibial insert,7#/8#, width 75mm,height 18mm	



NX-CR Insert

5218-10011	CR Tibial insert, 1#/2#, width 61mm,height 10mm	UHMWPE
5218-12011	CR Tibial insert, 1#/2#, width 61mm,height 12mm	
5218-14011	CR Tibial insert, 1#/2#, width 61mm,height 14mm	
5218-16011	CR Tibial insert, 1#/2#, width 61mm,height 16mm	
5218-18011	CR Tibial insert, 1#/2#, width 61mm,height 18mm	
5218-10021	CR Tibial insert, 3#/4#, width 66mm,height 10mm	
5218-12021	CR Tibial insert, 3#/4#, width 66mm,height 12mm	
5218-14021	CR Tibial insert, 3#/4#, width 66mm,height 14mm	
5218-16021	CR Tibial insert, 3#/4#, width 66mm,height 16mm	
5218-18021	CR Tibial insert, 3#/4#, width 66mm,height 18mm	
5218-10031	CR Tibial insert, 5#/6#, width 70mm,height 10mm	
5218-12031	CR Tibial insert, 5#/6#, width 70mm,height 12mm	
5218-14031	CR Tibial insert, 5#/6#, width 70mm,height 14mm	
5218-16031	CR Tibial insert, 5#/6#, width 70mm,height 16mm	
5218-18031	CR Tibial insert, 5#/6#, width 70mm,height 18mm	
5218-10041	CR Tibial insert, 7#/8#, width 75mm,height 10mm	
5218-12041	CR Tibial insert, 7#/8#, width 75mm,height 12mm	
5218-14041	CR Tibial insert, 7#/8#, width 75mm,height 14mm	
5218-16041	CR Tibial insert, 7#/8#, width 75mm,height 16mm	
5218-18041	CR Tibial insert, 7#/8#, width 75mm,height 18mm	



NX-Tibial Tray

Code	Description	Material
5216-01010	Tibial Tray, 1#, width 61mm	CoCrMo
5216-01020	Tibial Tray, 2#, width 63mm	
5216-01030	Tibial Tray, 3#, width 66mm	
5216-01040	Tibial Tray, 4#, width 68mm	
5216-01050	Tibial Tray, 5#, width 71mm	
5216-01060	Tibial Tray, 6#, width 73mm	
5216-01070	Tibial Tray, 7#, width 75mm	
5216-01080	Tibial Tray, 8#, width 77mm	

XM Primay Knee(PS&CR)





XM-PS Femoral Condyle

Code	Description	Material
5110-03021	PS Femoral Condyle (CLGMC) , 2#, left, width 57mm	CoCrMo
5110-03022	PS Femoral Condyle (CLGMC) , 2#, right, width 57mm	
5110-03031	PS Femoral Condyle (CLGMC) , 3#, left, width 60mm	
5110-03032	PS Femoral Condyle (CLGMC) , 3#, right, width 60mm	
5110-03041	PS Femoral Condyle (CLGMC) , 4#, left, width 62mm	
5110-03042	PS Femoral Condyle (CLGMC) , 4#, right, width 62mm	
5110-03051	PS Femoral Condyle (CLGMC) , 5#, left, width 65mm	
5110-03052	PS Femoral Condyle (CLGMC) , 5#, right, width 65mm	
5110-03061	PS Femoral Condyle (CLGMC) , 6#, left, width 67mm	
5110-03062	PS Femoral Condyle (CLGMC) , 6#, right, width 67mm	
5110-03071	PS Femoral Condyle (CLGMC) , 7#, left, width 70mm	
5110-03072	PS Femoral Condyle (CLGMC) , 7#, right, width 70mm	
5110-03081	PS Femoral Condyle (CLGMC) , 8#, left, width 72mm	
5110-03082	PS Femoral Condyle (CLGMC) , 8#, right, width 72mm	
5110-03101	PS Femoral Condyle (CLGMC) , 10#, left, width 77mm	
5110-03102	PS Femoral Condyle (CLGMC) , 10#, right, width 77mm	



XM-CR Femoral Condyle

Code	Description	Material
5114-03011	CR Femoral Condyle (CLGMC) , 2#, left, width 57mm	CoCrMo
5114-03012	CR Femoral Condyle (CLGMC) , 2#, right, width 57mm	
5114-03021	CR Femoral Condyle (CLGMC) , 3#, left, width 60mm	
5114-03022	CR Femoral Condyle (CLGMC) , 3#, right, width 60mm	
5114-03031	CR Femoral Condyle (CLGMC) , 4#, left, width 62mm	
5114-03032	CR Femoral Condyle (CLGMC) , 4#, right, width 62mm	
5114-03041	CR Femoral Condyle (CLGMC) , 5#, left, width 65mm	
5114-03042	CR Femoral Condyle (CLGMC) , 5#, right, width 65mm	
5114-03051	CR Femoral Condyle (CLGMC) , 6#, left, width 67mm	
5114-03052	CR Femoral Condyle (CLGMC) , 6#, right, width 67mm	
5114-03061	CR Femoral Condyle (CLGMC) , 7#, left, width 70mm	
5114-03062	CR Femoral Condyle (CLGMC) , 7#, right, width 70mm	
5114-03071	CR Femoral Condyle (CLGMC) , 8#, left, width 72mm	
5114-03072	CR Femoral Condyle (CLGMC) , 8#, right, width 72mm	
5114-03091	CR Femoral Condyle (CLGMC) , 10#, left, width 77mm	
5114-03092	CR Femoral Condyle (CLGMC) , 10#, right, width 77mm	



XM-PS Insert

Code	Description	Material
5112-09021	PS Tibial Insert (CLJMP) , 2#, width 62mm, height 9mm	UHMWPE
5112-11021	PS Tibial Insert (CLJMP) , 2#, width 62mm,height 11mm	
5112-13021	PS Tibial Insert (CLJMP) , 2#, width 62mm, height 13mm	
5112-15021	PS Tibial Insert (CLJMP) , 2#, width 62mm, height 15mm	
5112-17021	PS Tibial Insert (CLJMP) , 2#, width 62mm, height 17mm	
5112-19021	PS Tibial Insert (CLJMP) , 2#, width 62mm,height 19mm	
5112-09041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 9mm	
5112-11041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 11mm	
5112-13041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 13mm	
5112-15041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 15mm	
5112-17041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 17mm	
5112-19041	PS Tibial Insert (CLJMP) , 3/4#, width 67mm,height 19mm	
5112-09061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 9mm	
5112-11061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 11mm	
5112-13061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 13mm	
5112-15061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 15mm	
5112-17061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 17mm	
5112-19061	PS Tibial Insert (CLJMP) , 5/6#, width 72mm,height 19mm	
5112-09081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 9mm	
5112-11081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 11mm	

5112-13081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 13mm	
5112-15081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 15mm	
5112-17081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 17mm	
5112-19081	PS Tibial Insert (CLJMP) , 7/8#, width 77mm,height 19mm	
5112-09101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 9mm	
5112-11101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 11mm	
5112-13101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 13mm	
5112-15101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 15mm	
5112-17101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 17mm	
5112-19101	PS Tibial Insert (CLJMP) , 10#, width 82mm,height 19mm	

XM-CR Insert

Code	Description	Material
5115-09011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 9mm	UHMWPE
5115-11011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 11mm	
5115-13011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 13mm	
5115-15011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 15mm	
5115-17011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 17mm	
5115-19011	CR Tibial Insert (CLJMP) , 2#, width 62mm,height 19mm	
5115-09031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 9mm	
5115-11031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 11mm	
5115-13031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 13mm	
5115-15031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 15mm	
5115-17031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 17mm	
5115-19031	CR Tibial Insert (CLJMP) , 3/4#, width 67mm,height 19mm	
5115-09051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 9mm	
5115-11051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 11mm	
5115-13051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 13mm	
5115-15051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 15mm	

5115-17051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 17mm
5115-19051	CR Tibial Insert (CLJMP) , 5/6#, width 72mm,height 19mm
5115-09071	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 9mm
5115-11071	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 11mm
5115-13071	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 13mm
5115-15071	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 15mm
5115-17071	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 17mm
5115-19091	CR Tibial Insert (CLJMP) , 7/8#, width 77mm,height 19mm
5115-09091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 9mm
5115-11091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 11mm
5115-13091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 13mm
5115-15091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 15mm
5115-17091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 17mm
5115-19091	CR Tibial Insert (CLJMP) , 10#, width 82mm,height 19mm





XM-Tibial Tray

Code	Description	Material
5111-01020	Tbial Tray (CLJMC) , 2#, width 62mm	CoCrMo
5111-01030	Tibial Tray (CLJMC) , 3#, width 65mm	
5111-01040	Tibial Tray (CLJMC) , 4#, width 67mm	
5111-01050	Tibial Tray (CLJMC) , 5#, width 70mm	
5111-01060	Tibial Tray (CLJMC) , 6#, width 71mm	
5111-01070	Tibial Tray (CLJMC) , 7#, width 75mm	
5111-01080	Tibial Tray (CLJMC) , 8#, width 76mm	
5111-01100	Tibial Tray (CLJMC) , 10#, width 82mm	



Patella

Code	Description	Material
5213-08260	Patella, height 8mm, diameter 26mm	UHMWPE
5213-08280	Patella, height 8mm, diameter 29mm	
5213-08300	Patella, height 8mm, diameter 32mm	
5213-08320	Patella, height 8mm, diameter 35mm	
5213-10260	Patella, height 10mm, diameter 26mm	
5213-10280	Patella, height 10mm, diameter 29mm	
5213-10300	Patella, height 10mm, diameter 32mm	
5213-10320	Patella, height 10mm, diameter 35mm	
5213-12260	Patella, height 12mm, diameter 26mm	
5213-12280	Patella, height 12mm, diameter 29mm	
5213-12300	Patella, height 12mm, diameter 32mm	
5213-12320	Patella, height 12mm, diameter 35mm	

Hinged knee joint prosthesis



Femoral condyle

Code	Description	width	Material
5510-01021	Cemented femoral condyle (CLGMC) , 2#, L, width 60mm, Uncoated	60	CoCrMo
5510-01022	Cemented femoral condyle (CLGMC) , 2#, R, width 60mm, Uncoated	60	CoCrMo
5510-01031	Cemented femoral condyle (CLGMC) , 3#, L, width 62mm, Uncoated	62	CoCrMo
5510-01032	Cemented femoral condyle (CLGMC) , 3#, R, width 62mm, Uncoated	62	CoCrMo
5510-01041	Cemented femoral condyle (CLGMC) , 4#, L, width 65mm, Uncoated	65	CoCrMo
5510-01042	Cemented femoral condyle (CLGMC) , 4#, R, width 65mm, Uncoated	65	CoCrMo
5510-01051	Cemented femoral condyle (CLGMC) , 5#, L, width 68.5mm, Uncoated	68.5	CoCrMo
5510-01052	Cemented femoral condyle (CLGMC) , 5#, R, width 68.5mm, Uncoated	68.5	CoCrMo
5510-01061	Cemented femoral condyle (CLGMC) , 6#, L, width 71mm, Uncoated	71	CoCrMo
5510-01062	Cemented femoral condyle (CLGMC) , 6#, R, width 71mm, Uncoated	71	CoCrMo
5510-01071	Cemented femoral condyle (CLGMC) , 7#, L, width 75mm, Uncoated	75	CoCrMo
5510-01072	Cemented femoral condyle (CLGMC) , 7#, R, width 75mm, Uncoated	75	CoCrMo



Tibial tray

Code	Description	width	Material
5511-01020	Cemented tibial tray (CLJMC) , 2#, width 65mm, Uncoated	65	CoCrMo
5511-01030	Cemented tibial tray (CLJMC) , 3#, width 68mm, Uncoated	68	CoCrMo
5511-01040	Cemented tibial tray (CLJMC) , 4#, width 70mm, Uncoated	70	CoCrMo
5511-01050	Cemented tibial tray (CLJMC) , 5#, width 72mm, Uncoated	72	CoCrMo
5511-01060	Cemented tibial tray (CLJMC) , 6#, width 75mm, Uncoated	75	CoCrMo
5511-01070	Cemented tibial tray (CLJMC) , 7#, width 79mm, Uncoated	79	CoCrMo



Tibial insert

Code	Description	width	Material
5512-08021	Tibial insert (CLJMP) , 2#, thickness 8 mm, width 60mm	60	UHMWPE/XLPE
5512-10021	Tibial insert (CLJMP) , 2#, thickness 10mm, width 60mm	60	UHMWPE/XLPE
5512-12021	Tibial insert (CLJMP) , 2#, thickness 12mm, width 60mm	60	UHMWPE/XLPE
5512-14021	Tibial insert (CLJMP) , 2#, thickness 14 mm, width 60mm	60	UHMWPE/XLPE
5512-16021	Tibial insert (CLJMP) , 2#, thickness 16mm, width 60mm	60	UHMWPE/XLPE
5512-18021	Tibial insert (CLJMP) , 2#, thickness 18mm, width 60mm	60	UHMWPE/XLPE
5512-08031	Tibial insert (CLJMP) , 3#, thickness 8mm, width 62mm	62	UHMWPE/XLPE
5512-10031	Tibial insert (CLJMP) , 3#, thickness 10mm, width 62mm	62	UHMWPE/XLPE
5512-12031	Tibial insert (CLJMP) , 3#, thickness 12mm, width 62mm	62	UHMWPE/XLPE
5512-14031	Tibial insert (CLJMP) , 3#, thickness 14 mm, width 62mm	62	UHMWPE/XLPE
5512-16031	Tibial insert (CLJMP) , 3#, thickness 16mm, width 62mm	62	UHMWPE/XLPE
5512-18031	Tibial insert (CLJMP) , 3#, thickness 18mm, width 62mm	62	UHMWPE/XLPE
5512-08041	Tibial insert (CLJMP) , 4#, thickness 8mm, width 65mm	65	UHMWPE/XLPE
5512-10041	Tibial insert (CLJMP) , 4#, thickness 10mm, width 65mm	65	UHMWPE/XLPE
5512-12041	Tibial insert (CLJMP) , 4#, thickness 12mm, width 65mm	65	UHMWPE/XLPE
5512-14041	Tibial insert (CLJMP) , 4#, thickness 14 mm, width 65mm	65	UHMWPE/XLPE
5512-16041	Tibial insert (CLJMP) , 4#, thickness 16mm, width 65mm	65	UHMWPE/XLPE
5512-18041	Tibial insert (CLJMP) , 4#, thickness 18mm, width 65mm	65	UHMWPE/XLPE

5512-08051	Tibial insert (CLJMP) , 5#, thickness 8mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-10051	Tibial insert (CLJMP) , 5#, thickness 10mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-12051	Tibial insert (CLJMP) , 5#, thickness 12mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-14051	Tibial insert (CLJMP) , 5#, thickness 14 mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-16051	Tibial insert (CLJMP) , 5#, thickness 16mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-18051	Tibial insert (CLJMP) , 5#, thickness 18mm, width 68.5mm	68.5	UHMWPE/XLPE
5512-08061	Tibial insert (CLJMP) , 6#, thickness 8mm, width 71mm	71	UHMWPE/XLPE
5512-10061	Tibial insert (CLJMP) , 6#, thickness 10mm, width 71mm	71	UHMWPE/XLPE
5512-12061	Tibial insert (CLJMP) , 6#, thickness 12mm, width 71mm	71	UHMWPE/XLPE
5512-14061	Tibial insert (CLJMP) , 6#, thickness 14 mm, width 71mm	71	UHMWPE/XLPE
5512-16061	Tibial insert (CLJMP) , 6#, thickness 16mm, width 71mm	71	UHMWPE/XLPE
5512-18061	Tibial insert (CLJMP) , 6#, thickness 18mm, width 71mm	71	UHMWPE/XLPE
5512-08071	Tibial insert (CLJMP) , 7#, thickness 8mm, width 75mm	75	UHMWPE/XLPE
5512-10071	Tibial insert (CLJMP) , 7#, thickness 10mm, width 75mm	75	UHMWPE/XLPE
5512-12071	Tibial insert (CLJMP) , 7#, thickness 12mm, width 75mm	75	UHMWPE/XLPE
5512-14071	Tibial insert (CLJMP) , 7#, thickness 14 mm, width 75mm	75	UHMWPE/XLPE
5512-16071	Tibial insert (CLJMP) , 7#, thickness 16mm, width 75mm	75	UHMWPE/XLPE
5512-18071	Tibial insert (CLJMP) , 7#, thickness 18mm, width 75mm	75	UHMWPE/XLPE



Femoral extension

Code	Description	width	Material
5573-09030	Femoral extension, Diameter 9mm, Length 30mm	30	Titanium
5573-10030	Femoral extension, Diameter 10mm, Length 30mm	30	Titanium
5573-11030	Femoral extension, Diameter 11mm, Length 30mm	30	Titanium
5573-12030	Femoral extension, Diameter 12mm, Length 30mm	30	Titanium
5573-13030	Femoral extension, Diameter 13mm, Length 30mm	30	Titanium
5573-14030	Femoral extension, Diameter 14mm, Length 30mm	30	Titanium
5573-15030	Femoral extension, Diameter 15mm, Length 30mm	30	Titanium
5573-16030	Femoral extension, Diameter 16mm, Length 30mm	30	Titanium
5573-09060	Femoral extension, Diameter 9mm, Length 60mm	60	Titanium
5573-10060	Femoral extension, Diameter 10mm, Length 60mm	60	Titanium
5573-11060	Femoral extension, Diameter 11mm, Length 60mm	60	Titanium
5573-12060	Femoral extension, Diameter 12mm, Length 60mm	60	Titanium
5573-13060	Femoral extension, Diameter 13mm, Length 60mm	60	Titanium
5573-14060	Femoral extension, Diameter 14mm, Length 60mm	60	Titanium
5573-15060	Femoral extension, Diameter 15mm, Length 60mm	60	Titanium
5573-16060	Femoral extension, Diameter 16mm, Length 60mm	60	Titanium
5573-09100	Femoral extension, Diameter 9mm, Length 100mm	100	Titanium
5573-10100	Femoral extension, Diameter 10mm, Length 100mm	100	Titanium

5573-11100	Femoral extension, Diameter 11mm, Length 100mm	100	Titanium
5573-12100	Femoral extension, Diameter 12mm, Length 100mm	100	Titanium
5573-13100	Femoral extension, Diameter 13mm, Length 100mm	100	Titanium
5573-14100	Femoral extension, Diameter 14mm, Length 100mm	100	Titanium
5573-15100	Femoral extension, Diameter 15mm, Length 100mm	100	Titanium
5573-16100	Femoral extension, Diameter 16mm, Length 100mm	100	Titanium
5573-09063	Femoral extension with 3mm eccentricity, Diameter 9mm, Length 60mm	60	Titanium
5573-10063	Femoral extension with 3mm eccentricity, Diameter 10mm, Length 60mm	60	Titanium
5573-11063	Femoral extension with 3mm eccentricity, Diameter 11mm, Length 60mm	60	Titanium
5573-12063	Femoral extension with 3mm eccentricity, Diameter 12mm, Length 60mm	60	Titanium
5573-13083	Femoral extension with 3mm eccentricity, Diameter 13mm, Length 60mm	60	Titanium
5273-14083	Femoral extension with 3mm eccentricity, Diameter 14mm, Length 60mm	60	Titanium
5573-15083	Femoral extension with 3mm eccentricity, Diameter 15mm, Length 60mm	60	Titanium
5573-16083	Femoral extension with 3mm eccentricity, Diameter 16mm, Length 60mm	60	Titanium
5573-09103	Femoral extension with 3mm eccentricity, Diameter 9mm, Length 100mm	100	Titanium
5573-10103	Femoral extension with 3mm eccentricity, Diameter 10mm, Length 100mm	100	Titanium
5573-11103	Femoral extension with 3mm eccentricity, Diameter 11mm, Length 100mm	100	Titanium
5573-12103	Femoral extension with 3mm eccentricity, Diameter 12mm, Length 100mm	100	Titanium
5573-13103	Femoral extension with 3mm eccentricity, Diameter 13mm, Length 100mm	100	Titanium
5273-14103	Femoral extension with 3mm eccentricity, Diameter 14mm, Length 100mm	100	Titanium
5573-15103	Femoral extension with 3mm eccentricity, Diameter 15mm, Length 100mm	100	Titanium
5573-16103	Femoral extension with 3mm eccentricity, Diameter 16mm, Length 100mm	100	Titanium



Tibial extension

Code	Description	width	Material
5574-09030	Tibial extension, Diameter 9mm, Length 30mm	30	Titanium
5574-10030	Tibial extension, Diameter 10mm, Length 30mm	30	Titanium
5574-11030	Tibial extension, Diameter 11mm, Length 30mm	30	Titanium
5574-12030	Tibial extension, Diameter 12mm, Length 30mm	30	Titanium
5574-13030	Tibial extension, Diameter 13mm, Length 30mm	30	Titanium
5574-14030	Tibial extension, Diameter 14mm, Length 30mm	30	Titanium
5574-15030	Tibial extension, Diameter 15mm, Length 30mm	30	Titanium
5574-16030	Tibial extension, Diameter 16mm, Length 30mm	30	Titanium
5574-09060	Tibial extension, Diameter 9mm, Length 60mm	60	Titanium
5574-10060	Tibial extension, Diameter 10mm, Length 60mm	60	Titanium
5574-11060	Tibial extension, Diameter 11mm, Length 60mm	60	Titanium
5574-12060	Tibial extension, Diameter 12mm, Length 60mm	60	Titanium
5574-13060	Tibial extension, Diameter 13mm, Length 60mm	60	Titanium
5574-14060	Tibial extension, Diameter 14mm, Length 60mm	60	Titanium
5574-15060	Tibial extension, Diameter 15mm, Length 60mm	60	Titanium
5574-16060	Tibial extension, Diameter 16mm, Length 60mm	60	Titanium
5574-09100	Tibial extension, Diameter 9mm, Length 100mm	100	Titanium
5574-10100	Femoral extension, Diameter 10mm, Length 100mm	100	Titanium
5574-11100	Tibial extension, Diameter 11mm, Length 100mm	100	Titanium

5574-12100	Tibial extension, Diameter 12mm, Length 100mm	100	Titanium
5574-13100	Tibial extension, Diameter 13mm, Length 100mm	100	Titanium
5574-14100	Tibial extension, Diameter 14mm, Length 100mm	100	Titanium
5574-15100	Tibial extension, Diameter 15mm, Length 100mm	100	Titanium
5574-16100	Tibial extension, Diameter 16mm, Length 100mm	100	Titanium
5574-09063	Tibial extension with 3mm eccentricity, Diameter 9mm, Length 60mm	60	Titanium
5574-10063	Tibial extension with 3mm eccentricity, Diameter 10mm, Length 60mm	60	Titanium
5574-11063	Tibial extension with 3mm eccentricity, Diameter 11mm, Length 60mm	60	Titanium
5574-12063	Tibial extension with 3mm eccentricity, Diameter 12mm, Length 60mm	60	Titanium
5574-13083	Tibial extension with 3mm eccentricity, Diameter 13mm, Length 60mm	60	Titanium
5274-14083	Tibial extension with 3mm eccentricity, Diameter 14mm, Length 60mm	60	Titanium
5574-15083	Tibial extension with 3mm eccentricity, Diameter 15mm, Length 60mm	60	Titanium
5574-16083	Tibial extension with 3mm eccentricity, Diameter 16mm, Length 60mm	60	Titanium
5574-09103	Tibial extension with 3mm eccentricity, Diameter 9mm, Length 100mm	100	Titanium
5574-10103	Tibial extension with 3mm eccentricity, Diameter 10mm, Length 100mm	100	Titanium
5574-11103	Tibial extension with 3mm eccentricity, Diameter 11mm, Length 100mm	100	Titanium
5574-12103	Tibial extension with 3mm eccentricity, Diameter 12mm, Length 100mm	100	Titanium
5574-13103	Tibial extension with 3mm eccentricity, Diameter 13mm, Length 100mm	100	Titanium
5274-14103	Tibial extension with 3mm eccentricity, Diameter 14mm, Length 100mm	100	Titanium
5574-15103	Tibial extension with 3mm eccentricity, Diameter 15mm, Length 100mm	100	Titanium
5574-16103	Tibial extension with 3mm eccentricity, Diameter 16mm, Length 100mm	100	Titanium

Wedge

Femoral wedge



Femoral Condyle- Distal wedge

Code	Description	width	Material
5572-05023	2#-5mm	5mm	Titanium
5572-10023	2#-10mm	10mm	Titanium
5572-05033	3#-5mm	5mm	Titanium
5572-10033	3#-10mm	10mm	Titanium
5572-05043	4#-5mm	5mm	Titanium
5572-10043	4#-10mm	10mm	Titanium
5572-05053	5#-5mm	5mm	Titanium
5572-10053	5#-10mm	10mm	Titanium
5572-05063	6#-5mm	5mm	Titanium
5572-10063	6#-10mm	10mm	Titanium
5572-05073	7#-5mm	5mm	Titanium
5572-10073	7#-10mm	10mm	Titanium



Femoral Condyle- Posterior wedge

Code	Description	width	Material
5571-05023	2#-5mm	5mm	Titanium
5571-10023	2#-10mm	10mm	Titanium
5571-05033	3#-5mm	5mm	Titanium
5571-10033	3#-10mm	10mm	Titanium
5571-05043	4#-5mm	5mm	Titanium
5571-10043	4#-10mm	10mm	Titanium
5571-05053	5#-5mm	5mm	Titanium
5571-10053	5#-10mm	10mm	Titanium
5571-05063	6#-5mm	5mm	Titanium
5571-10063	6#-10mm	10mm	Titanium
5571-05073	7#-5mm	5mm	Titanium
5571-10073	7#-10mm	10mm	Titanium

Tibial wedge



Tibial wedge (Half-length wedge)

Code	Description	width	Material
5576-05021	2#L-5mm	5mm	Titanium
5576-10021	2#L-10mm	10mm	Titanium
5576-05031	3#L-5mm	5mm	Titanium
5576-10031	3#L-10mm	10mm	Titanium
5576-05041	4#L-5mm	5mm	Titanium
5576-10041	4#L-10mm	10mm	Titanium
5576-05051	5#L-5mm	5mm	Titanium
5576-10051	5#L-10mm	10mm	Titanium
5576-05061	6#L-5mm	5mm	Titanium
5576-10061	6#L-10mm	10mm	Titanium
5576-05071	7#L-5mm	5mm	Titanium
5576-10071	7#L-10mm	10mm	Titanium
5576-05022	2#R-5mm	5mm	Titanium
5576-10022	2#R-10mm	10mm	Titanium
5576-05032	3#R-5mm	5mm	Titanium
5576-10032	3#R-10mm	10mm	Titanium
5576-05042	4#R-5mm	5mm	Titanium

5576-10042	4#R-10mm	10mm	Titanium
5576-05052	5#R-5mm	5mm	Titanium
5576-10052	5#R-10mm	10mm	Titanium
5576-05062	6#R-5mm	5mm	Titanium
5576-10062	6#R-10mm	10mm	Titanium
5576-05072	7#R-5mm	5mm	Titanium
5576-10072	7#R-10mm	10mm	Titanium



Tibial wedge (Half-length cushion wedge)

Code	Description	width	Material
5577-10021	2#L-10mm	10mm	Titanium
5577-10031	3#L-10mm	10mm	Titanium
5577-10041	4#L-10mm	10mm	Titanium
5577-10051	5#L-10mm	10mm	Titanium
5577-10061	6#L-10mm	10mm	Titanium
5577-10071	7#L-10mm	10mm	Titanium
5577-10022	2#R-10mm	10mm	Titanium
5577-10032	3#R-10mm	10mm	Titanium
5577-10042	4#R-10mm	10mm	Titanium
5577-10052	5#R-10mm	10mm	Titanium
5577-10062	6#R-10mm	10mm	Titanium
5577-10072	7#R-10mm	10mm	Titanium



Tibial wedge (Full-length cushion wedge)

Code	Description	width	Material
5578-10021	2#L-10mm	10mm	Titanium
5578-10031	3#L-10mm	10mm	Titanium
5578-10041	4#L-10mm	10mm	Titanium
5578-10051	5#L-10mm	10mm	Titanium
5578-10061	6#L-10mm	10mm	Titanium
5578-10071	7#L-10mm	10mm	Titanium
5578-10022	2#R-10mm	10mm	Titanium
5578-10032	3#R-10mm	10mm	Titanium
5578-10042	4#R-10mm	10mm	Titanium
5578-10052	5#R-10mm	10mm	Titanium
5578-10062	6#R-10mm	10mm	Titanium
5578-10072	7#R-10mm	10mm	Titanium



Wing sleeve

Code	Description	width	Material
5579-28000	28mm	28mm	Titanium

NX/CCK System



NX Insert



NX Insert




NX Insert





NX-Femoral Condyle

Code	Description	Material
5270-01021	Femoral Condyle (CLGMC) 2#, left, width 57mm	CoCrMo
5270-01022	Femoral Condyle (CLGMC) 2#, right, width 57mm	
5270-01031	Femoral Condyle (CLGMC) 3#, left, width 60mm	
5270-01032	Femoral Condyle (CLGMC) 3#, right, width 60mm	
5270-01041	Femoral Condyle (CLGMC) 4#, left, width 62mm	
5270-01042	Femoral Condyle (CLGMC) 4#, right, width 62mm	
5270-01051	Femoral Condyle (CLGMC) 5#, left, width 65mm	
5270-01052	Femoral Condyle (CLGMC) 5#, right, width 65mm	
5270-01061	Femoral Condyle (CLGMC) 6#, left, width 67mm	
5270-01062	Femoral Condyle (CLGMC) 6#, right, width 67mm	
5270-01071	Femoral Condyle (CLGMC) 7#, left, width 70mm	
5270-01072	Femoral Condyle (CLGMC) 7#, right, width 70mm	



NX-Insert

Code	Description	Material
5271-09031	Tibial Insert (CLJMP),3#/4#,Height 9mm,width 68mm	UHMWPE
5271-11031	Tibial Insert (CLJMP),3#/4#,Height 11mm,width 68mm	
5271-13031	Tibial Insert (CLJMP),3#/4#,Height 13mm,width 68mm	
5271-15031	Tibial Insert (CLJMP),3#/4#,Height 15mm,width 68mm	
5271-17031	Tibial Insert (CLJMP),3#/4#,Height 17mm,width 68mm	
5271-19031	Tibial Insert (CLJMP),3#/4#,Height 19mm,width 68mm	
5271-09051	Tibial Insert (CLJMP),5#/6#,Height 9mm,width 73mm	
5271-11051	Tibial Insert (CLJMP),5#/6#,Height 11mm,width 73mm	
5271-13051	Tibial Insert (CLJMP),5#/6#,Height 13mm,width 73mm	
5271-15051	Tibial Insert (CLJMP),5#/6#,Height 15mm,width 73mm	
5271-17051	Tibial Insert (CLJMP),5#/6#,Height 17mm,width 73mm	
5271-19051	Tibial Insert (CLJMP),5#/6#,Height 19mm,width 73mm	
5271-09071	Tibial Insert (CLJMP),7#,Height 9mm,width 78mm	
5271-11071	Tibial Insert (CLJMP),7#,Height 11mm,width 78mm	
5271-13071	Tibial Insert (CLJMP),7#,Height 13mm,width 78mm	
5271-15071	Tibial Insert (CLJMP),7#,Height 15mm,width 78mm	
5271-17071	Tibial Insert (CLJMP),7#,Height 17mm,width 78mm	
5271-19071	Tibial Insert (CLJMP),7#,Height 19mm,width 78mm	



NX- Tibial Tray

Code	Description	Material
5171-01020	Tibial tray(CLJMC) 2#, width 61mm	CoCrMo
5171-01030	Tibial tray(CLJMC) 3#, width 65mm	
5171-01040	Tibial tray(CLJMC) 4#, width 68mm	
5171-01050	Tibial tray(CLJMC) 5#, width 70mm	
5171-01060	Tibial tray(CLJMC) 6#, width 73mm	
5171-01070	Tibial tray(CLJMC) 7#, width 78mm	



Straight Femoral Stem

Code	Description	Material
/	/	Ti6Al4V
5273-10060	Femoral Stem, diameter10mm, length 80mm	
5273-11060	Femoral Stem, diameter11mm, length 80mm	
5273-12060	Femoral Stem, diameter12mm, length 80mm	
5273-13060	Femoral Stem, diameter13mm, length 80mm	
5273-14060	Femoral Stem, diameter14mm, length 80mm	
5273-15060	Femoral Stem, diameter15mm, length 80mm	
/	/	
/	/	
5273-10100	Femoral Stem, diameter10mm, length 100mm	
5273-11100	Femoral Stem, diameter11mm, length 100mm	
5273-12100	Femoral Stem, diameter12mm, length 100mm	
5273-13100	Femoral Stem, diameter13mm, length 100mm	
5273-14100	Femoral Stem, diameter14mm, length 100mm	
5273-15100	Femoral Stem, diameter15mm, length 100mm	
5273-16100	Femoral Stem, diameter16mm, length 100mm	
/	/	
5273-10130	Femoral Stem, diameter10mm, length 120mm	
5273-11130	Femoral Stem, diameter11mm, length 120mm	
5273-12130	Femoral Stem, diameter12mm, length 120mm	
5273-13130	Femoral Stem, diameter13mm, length 120mm	
5273-14130	Femoral Stem, diameter14mm, length 120mm	
5273-15130	Femoral Stem, diameter15mm, length 120mm	
/	/	



Femoral Stem with Offset

Code	Description	Material
/	/	Ti6Al4V
5273-10103	Femoral Stem with offset, diameter10mm, length 100mm	
5273-11103	Femoral Stem with offset, diameter11mm, length 100mm	
5273-12103	Femoral Stem with offset, diameter12mm, length 100mm	
5273-13103	Femoral Stem with offset, diameter13mm, length 100mm	
5273-14103	Femoral Stem with offset, diameter14mm, length 100mm	
5273-15103	Femoral Stem with offset, diameter15mm, length 100mm	
5273-16103	Femoral Stem with offset, diameter16mm, length 100mm	
/	/	
5273-10133	Femoral Stem with offset, diameter10mm, length 120mm	
5273-11133	Femoral Stem with offset, diameter11mm, length 120mm	
5273-12133	Femoral Stem with offset, diameter12mm, length 120mm	
5273-13133	Femoral Stem with offset, diameter13mm, length 120mm	
5273-14133	Femoral Stem with offset, diameter14mm, length 120mm	
5273-15133	Femoral Stem with offset, diameter15mm, length 120mm	
/	/	



Fixation Screw

Code	Description	Material
5275-10050	Fixation Screw,Diameter 5mm	Ti6Al4V
5275-10060	Fixation Screw,Diameter 6mm	
5275-10080	Fixation Screw,Diameter 8mm	

Knee prosthesis



Knee Joint Prosthesis-New Axial

Code	Description	width	Material
6240-01021	Femoral New Axis 2#, L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6240-01022	Femoral New Axis 2#, R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6240-01041	Femoral New Axis 4#, L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6240-01042	Femoral New Axis 4#, R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6240-01061	Femoral New Axis 6#, L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6240-01062	Femoral New Axis 6#, R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6240-01081	Femoral New Axis 8#, L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6240-01082	Femoral New Axis 8#, R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6240-01101	Femoral New Axis 10#, L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6240-01102	Femoral New Axis 10#, R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6241-01021	Tibial New Axis 2#, L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6241-01022	Tibial New Axis 2#, R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6241-01041	Tibial New Axis 4#, L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6241-01042	Tibial New Axis 4#, R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6241-01061	Tibial New Axis 6#, L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE

6241-01062	Tibial New Axis 6#, R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6241-01081	Tibial New Axis 8#, L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6241-01082	Tibial New Axis 8#, R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6241-01101	Tibial New Axis 10#, L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6241-01102	Tibial New Axis 10#, R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE



Knee prosthesis-stem

6270-45131	Femoral special stem13/145/10, Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium
6270-45231	Femoral special stem13/145/20Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium
6270-55131	Femoral special stem13/155/10Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium
6270-55231	Femoral special stem13/155/20Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium
6270-65131	Femoral special stem13/165/10Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium
6270-65231	Femoral special stem13/165/20Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium
6270-45141	Femoral special stem14/145/10Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium
6270-45241	Femoral special stem14/145/20Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium

6270-55141	Femoral special stem14/155/10Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium,
6270-55241	Femoral special stem14/155/20Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium,
6270-65141	Femoral special stem14/165/10Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium,
6270-65241	Femoral special stem14/165/20Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium,
6270-45151	Femoral special stem15/145/10Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium,
6270-45251	Femoral special stem15/145/20Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium,
6270-55151	Femoral special stem15/155/10Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium,
6270-55251	Femoral special stem15/155/20Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium,
6270-65151	Femoral special stem15/165/10Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium,
6270-65251	Femoral special stem15/165/20Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium,
6271-35311	Universal Stem 11/135/30Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35411	Universal Stem 11/135/40Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35511	Universal Stem 11/135/50Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,

6271-35611	Universal Stem 11/135/60Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45311	Universal Stem 11/145/30Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45411	Universal Stem 11/145/40Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45511	Universal Stem 11/145/50Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45611	Universal Stem 11/145/60Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55311	Universal Stem 11/155/30Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55411	Universal Stem 11/155/40Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55511	Universal Stem 11/155/50Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55611	Universal Stem 11/155/60Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65311	Universal Stem 11/165/30Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65411	Universal Stem 11/165/40Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65511	Universal Stem 11/165/50Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65611	Universal Stem 11/165/60Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,

6271-35321	Universal Stem 12/135/30Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35421	Universal Stem 12/135/40Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35521	Universal Stem 12/135/50Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35621	Universal Stem 12/135/60Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45321	Universal Stem 12/145/30Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45421	Universal Stem 12/145/40Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45521	Universal Stem 12/145/50Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45621	Universal Stem 12/145/60Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55321	Universal Stem 12/155/30Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55421	Universal Stem 12/155/40Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55521	Universal Stem 12/155/50Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55621	Universal Stem 12/155/60Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65321	Universal Stem 12/165/30Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,

6271-65421	Universal Stem 12/165/40Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65521	Universal Stem 12/165/50Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65621	Universal Stem 12/165/50Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35331	Universal Stem 13/135/30Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35431	Universal Stem 13/135/40Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35531	Universal Stem 13/135/50Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35631	Universal Stem 13/135/60Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45331	Universal Stem 13/145/30Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45431	Universal Stem 13/145/40Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45531	Universal Stem 13/145/50Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45631	Universal Stem 13/145/60Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55331	Universal Stem 13/155/30Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55431	Universal Stem 13/155/40Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,

6271-55531	Universal Stem 13/155/50Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55631	Universal Stem 13/155/60Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65331	Universal Stem 13/165/30Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65431	Universal Stem 13/165/40Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65531	Universal Stem 13/165/50Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65631	Universal Stem 13/165/60Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35341	Universal Stem 14/135/30Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35441	Universal Stem 14/135/40Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35541	Universal Stem 14/135/50Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35641	Universal Stem 14/135/60Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45341	Universal Stem 14/145/30Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45441	Universal Stem 14/145/40Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45541	Universal Stem 14/145/50Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,

6271-45641	Universal Stem 14/145/60Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55341	Universal Stem 14/155/30Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55441	Universal Stem 14/155/40Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55541	Universal Stem 14/155/50Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55641	Universal Stem 14/155/60Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65341	Universal Stem 14/165/30Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65441	Universal Stem 14/165/40Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65541	Universal Stem 14/165/50Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65641	Universal Stem 14/165/60Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35351	Universal Stem 15/135/30Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35451	Universal Stem 15/135/40Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35551	Universal Stem 15/135/50Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35651	Universal Stem 15/135/60Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,

6271-45351	Universal Stem 15/145/30Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45451	Universal Stem 15/145/40Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45551	Universal Stem 15/145/50Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45651	Universal Stem 15/145/60Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55351	Universal Stem 15/155/30Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55451	Universal Stem 15/155/40Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55551	Universal Stem 15/155/50Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55651	Universal Stem 15/155/60Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65351	Universal Stem 15/165/30Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65451	Universal Stem 15/165/40Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65551	Universal Stem 15/165/50Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65651	Universal Stem 15/165/60Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,



Knee prosthesis-Osteotomy segment

6272-24101	Osteotomy segment 24/100, Diameter 24 mm, Length 100 mm	100	Titanium,
6272-24111	Osteotomy segment 24/110, Diameter 24 mm, Length 110 mm	110	Titanium,
6272-24121	Osteotomy segment 24/120, Diameter 24 mm, Length 120 mm	120	Titanium,
6272-24041	Osteotomy segment 24/40, Diameter 24 mm, Length 40 mm	40	Titanium,
6272-24051	Osteotomy segment 24/50, Diameter 24 mm, Length 50 mm	50	Titanium,
6272-24061	Osteotomy segment 24/60, Diameter 24 mm, Length 60 mm	60	Titanium,
6272-24071	Osteotomy segment 24/70, Diameter 24 mm, Length 70 mm	70	Titanium,
6272-24081	Osteotomy segment 24/80, Diameter 24 mm, Length 80 mm	80	Titanium,
6272-24091	Osteotomy segment 24/90, Diameter 24 mm, Length 90 mm	90	Titanium,



Knee joint prosthesis-new axis (combined with ordinary)

Code	Description	width	Material
6240-02021	Femoral New Axis 2#, L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6240-02022	Femoral New Axis 2#, R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6240-02041	Femoral New Axis 4#, L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6240-02042	Femoral New Axis 4#, R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6240-02061	Femoral New Axis 6#, L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6240-02062	Femoral New Axis 6#, R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6240-02081	Femoral New Axis 8#, L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6240-02082	Femoral New Axis 8#, R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6240-02101	Femoral New Axis 10#, L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6240-02102	Femoral New Axis 10#, R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6241-02021	Tibial New Axis 2#, L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6241-02022	Tibial New Axis 2#, R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6241-02041	Tibial New Axis 4#, L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6241-02042	Tibial New Axis 4#, R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE

6241-02061	Tibial New Axis 6#, L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6241-02062	Tibial New Axis 6#, R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6241-02081	Tibial New Axis 8#, L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6241-02082	Tibial New Axis 8#, R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6241-02101	Tibial New Axis 10#, L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6241-02102	Tibial New Axis 10#, R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
Knee joint prosthesis-stem			
6270-45132	Femoral special stem13/145/10, Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium
6270-45232	Femoral special stem13/145/20Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium
6270-55132	Femoral special stem13/155/10Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium
6270-55232	Femoral special stem13/155/20Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium
6270-65132	Femoral special stem13/165/10Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium,
6270-65232	Femoral special stem13/165/20Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium,
6270-45142	Femoral special stem14/145/10Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium,
6270-45242	Femoral special stem14/145/20Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium,
6270-55142	Femoral special stem14/155/10Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium,

6270-55242	Femoral special stem14/155/20Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium,
6270-65142	Femoral special stem14/165/10Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium,
6270-65242	Femoral special stem14/165/20Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium,
6270-45152	Femoral special stem15/145/10Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 10 mm	145	Titanium,
6270-45252	Femoral special stem15/145/20Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 20 mm	145	Titanium,
6270-55152	Femoral special stem15/155/10Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 10 mm	155	Titanium,
6270-55252	Femoral special stem15/155/20Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 20 mm	155	Titanium,
6270-65152	Femoral special stem15/165/10Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 10 mm	165	Titanium,
6270-65252	Femoral special stem15/165/20Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 20 mm	165	Titanium,
6271-35312	Universal Stem 11/135/30Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35412	Universal Stem 11/135/40Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35512	Universal Stem 11/135/50Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35612	Universal Stem 11/135/60Stem Diameter 11 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,

6271-45312	Universal Stem 11/145/30Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45412	Universal Stem 11/145/40Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45512	Universal Stem 11/145/50Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45612	Universal Stem 11/145/60Stem Diameter 11 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55312	Universal Stem 11/155/30Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55412	Universal Stem 11/155/40Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55512	Universal Stem 11/155/50Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55612	Universal Stem 11/155/60Stem Diameter 11 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65312	Universal Stem 11/165/30Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65412	Universal Stem 11/165/40Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65512	Universal Stem 11/165/50Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65612	Universal Stem 11/165/60Stem Diameter 11 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35322	Universal Stem 12/135/30Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,

6271-35422	Universal Stem 12/135/40Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35522	Universal Stem 12/135/50Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35622	Universal Stem 12/135/60Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45322	Universal Stem 12/145/30Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45422	Universal Stem 12/145/40Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45522	Universal Stem 12/145/50Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45622	Universal Stem 12/145/60Stem Diameter 12 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55322	Universal Stem 12/155/30Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55422	Universal Stem 12/155/40Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55522	Universal Stem 12/155/50Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55622	Universal Stem 12/155/60Stem Diameter 12 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65322	Universal Stem 12/165/30Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65422	Universal Stem 12/165/40Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,

6271-65522	Universal Stem 12/165/50Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65622	Universal Stem 12/165/50Stem Diameter 12 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35332	Universal Stem 13/135/30Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35432	Universal Stem 13/135/40Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35532	Universal Stem 13/135/50Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35632	Universal Stem 13/135/60Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45332	Universal Stem 13/145/30Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45432	Universal Stem 13/145/40Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45532	Universal Stem 13/145/50Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45632	Universal Stem 13/145/60Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55332	Universal Stem 13/155/30Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55432	Universal Stem 13/155/40Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55532	Universal Stem 13/155/50Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,

6271-55632	Universal Stem 13/155/60Stem Diameter 13 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65332	Universal Stem 13/165/30Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65432	Universal Stem 13/165/40Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65532	Universal Stem 13/165/50Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65632	Universal Stem 13/165/60Stem Diameter 13 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35342	Universal Stem 14/135/30Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35442	Universal Stem 14/135/40Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35542	Universal Stem 14/135/50Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35642	Universal Stem 14/135/60Stem Diameter 14 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45342	Universal Stem 14/145/30Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45442	Universal Stem 14/145/40Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45542	Universal Stem 14/145/50Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45642	Universal Stem 14/145/60Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,

6271-55342	Universal Stem 14/155/30Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55442	Universal Stem 14/155/40Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55542	Universal Stem 14/155/50Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55642	Universal Stem 14/155/60Stem Diameter 14 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65342	Universal Stem 14/165/30Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65442	Universal Stem 14/165/40Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65542	Universal Stem 14/165/50Stem Diameter 14 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65642	Universal Stem 14/165/60Stem Diameter 14 mm, Stem Length 145 mm, With osteotomy segment 60 mm	165	Titanium,
6271-35352	Universal Stem 15/135/30Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35452	Universal Stem 15/135/40Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35552	Universal Stem 15/135/50Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35652	Universal Stem 15/135/60Stem Diameter 15 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45352	Universal Stem 15/145/30Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,

6271-45452	Universal Stem 15/145/40Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 40 mm	145	Titanium,
6271-45552	Universal Stem 15/145/50Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45652	Universal Stem 15/145/60Stem Diameter 15 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,
6271-55352	Universal Stem 15/155/30Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 30 mm	155	Titanium,
6271-55452	Universal Stem 15/155/40Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 40 mm	155	Titanium,
6271-55552	Universal Stem 15/155/50Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 50 mm	155	Titanium,
6271-55652	Universal Stem 15/155/60Stem Diameter 15 mm, Stem Length 155 mm, With osteotomy segment 60 mm	155	Titanium,
6271-65352	Universal Stem 15/165/30Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 30 mm	165	Titanium,
6271-65452	Universal Stem 15/165/40Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 40 mm	165	Titanium,
6271-65552	Universal Stem 15/165/50Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 50 mm	165	Titanium,
6271-65652	Universal Stem 15/165/60Stem Diameter 15 mm, Stem Length 165 mm, With osteotomy segment 60 mm	165	Titanium,
Knee joint prosthesis-Osteotomy segment			
6272-24102	Osteotomy segment 24/100, Diameter 24 mm, Length 100 mm	100	Titanium,
6272-24112	Osteotomy segment 24/110, Diameter 24 mm, Length 110 mm	110	Titanium,
6272-24122	Osteotomy segment 24/120, Diameter 24 mm, Length 120 mm	120	Titanium,

6272-24042	Osteotomy segment 24/40, Diameter 24 mm, Length 40 mm	40	Titanium,
6272-24052	Osteotomy segment 24/50, Diameter 24 mm, Length 50 mm	50	Titanium,
6272-24062	Osteotomy segment 24/60, Diameter 24 mm, Length 60 mm	60	Titanium,
6272-24072	Osteotomy segment 24/70, Diameter 24 mm, Length 70 mm	70	Titanium
6272-24082	Osteotomy segment 24/80, Diameter 24 mm, Length 80 mm	80	Titanium
6272-24092	Osteotomy segment 24/90, Diameter 24 mm, Length 90 mm	90	Titanium



Knee Joint Prosthesis-New Axial (Customized Aseptic)

Code	Description	width	Material
6640-01021	Femoral Adjustable New Axis 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6640-01041	Femoral Adjustable New Axis 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6640-01061	Femoral Adjustable New Axis 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6640-01081	Femoral Adjustable New Axis 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6640-01101	Femoral Adjustable New Axis 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6640-01022	Femoral Adjustable New Axis 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6640-01042	Femoral Adjustable New Axis 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6640-01062	Femoral Adjustable New Axis 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6640-01082	Femoral Adjustable New Axis 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6640-01102	Femoral Adjustable New Axis 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6642-01021	Tibial Adjustable New Axis 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6642-01041	Tibial Adjustable New Axis 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6642-01061	Tibial Adjustable New Axis 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6642-01081	Tibial Adjustable New Axis 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE

6642-01101	Tibial Adjustable New Axis 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6642-01022	Tibial Adjustable New Axis 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6642-01042	Tibial Adjustable New Axis 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6642-01062	Tibial Adjustable New Axis 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6642-01082	Tibial Adjustable New Axis 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6642-01102	Tibial Adjustable New Axis 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6644-01021	Child Reserved Epiphyseal Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6644-01041	Child Reserved Epiphyseal Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6644-01061	Child Reserved Epiphyseal Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6644-01081	Child Reserved Epiphyseal Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6644-01101	Child Reserved Epiphyseal Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6644-01022	Child Reserved Epiphyseal Axis, 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6644-01042	Child Reserved Epiphyseal Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6644-01062	Child Reserved Epiphyseal Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6644-01082	Child Reserved Epiphyseal Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6644-01102	Child Reserved Epiphyseal Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6645-01021	Femoral Adjustable Half Joint, 2L, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6645-01041	Femoral Adjustable Half Joint, 4L, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6645-01061	Femoral Adjustable Half Joint, 6L, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6645-01081	Femoral Adjustable Half Joint, 8L, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6645-01101	Femoral Adjustable Half Joint, 10L, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6645-01022	Femoral Adjustable Half Joint, 2R, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6645-01042	Femoral Adjustable Half Joint, 4R, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6645-01062	Femoral Adjustable Half Joint, 6R, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6645-01082	Femoral Adjustable Half Joint, 8R, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE

6645-01102	Femoral Adjustable Half Joint, 10R, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6646-01021	Tibial Adjustable Half Joint, 2LTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6646-01041	Tibial Adjustable Half Joint, 4LTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6646-01061	Tibial Adjustable Half Joint, 6LTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6646-01081	Tibial Adjustable Half Joint, 8LTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6646-01101	Tibial Adjustable Half Joint, 10LTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6646-01022	Tibial Adjustable Half Joint, 2RTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6646-01042	Tibial Adjustable Half Joint, 4RTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6646-01062	Tibial Adjustable Half Joint, 6RTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6646-01082	Tibial Adjustable Half Joint, 8RTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6646-01102	Tibial Adjustable Half Joint, 10RTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6647-01021	Tibia, Total Femoral Adjustable New Axis, 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6647-01041	Tibia, Total Femoral Adjustable New Axis, 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6647-01061	Tibia, Total Femoral Adjustable New Axis, 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6647-01081	Tibia, Total Femoral Adjustable New Axis, 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6647-01101	Tibia, Total Femoral Adjustable New Axis, 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6647-01022	Tibia, Total Femoral Adjustable New Axis, 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6647-01042	Tibia, Total Femoral Adjustable New Axis, 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6647-01062	Tibia, Total Femoral Adjustable New Axis, 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6647-01082	Tibia, Total Femoral Adjustable New Axis, 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6647-01102	Tibia, Total Femoral Adjustable New Axis, 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6648-01021	Femoral Epiphysis Reserved Adjustable Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6648-01041	Femoral Epiphysis Reserved Adjustable Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6648-01061	Femoral Epiphysis Reserved Adjustable Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6648-01081	Femoral Epiphysis Reserved Adjustable Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE

6648-01101	Femoral Epiphysis Reserved Adjustable Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6648-01022	Femoral Epiphysis Reserved Adjustable Axis, 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6648-01042	Femoral Epiphysis Reserved Adjustable Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6648-01062	Femoral Epiphysis Reserved Adjustable Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6648-01082	Femoral Epiphysis Reserved Adjustable Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6648-01102	Femoral Epiphysis Reserved Adjustable Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6658-01021	Half Femoral Joint, 2L, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6658-01041	Half Femoral Joint, 4L, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6658-01061	Half Femoral Joint, 6L, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6658-01081	Half Femoral Joint, 8L, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6658-01101	Half Femoral Joint, 10L, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6658-01022	Half Femoral Joint, 2R, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6658-01042	Half Femoral Joint, 4R, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6658-01062	Half Femoral Joint, 6R, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6658-01082	Half Femoral Joint, 8R, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6658-01102	Half Femoral Joint, 10R, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6659-01021	Half Tibial Joint, 2LTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6659-01041	Half Tibial Joint, 4LTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6659-01061	Half Tibial Joint, 6LTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6659-01081	Half Tibial Joint, 8LTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6659-01101	Half Tibial Joint, 10LTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6659-01022	Half Tibial Joint, 2RTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6659-01042	Half Tibial Joint, 4RTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6659-01062	Half Tibial Joint, 6RTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6659-01082	Half Tibial Joint, 8RTray width 72 mm	72	Titanium, CoCrMo, UHMWPE

6659-01102	Half Tibial Joint, 10R Tray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6660-01021	Femoral, Tibial New Axis , 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6660-01041	Femoral, Tibial New Axis , 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6660-01061	Femoral, Tibial New Axis , 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6660-01081	Femoral, Tibial New Axis , 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6660-01101	Femoral, Tibial New Axis , 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6660-01022	Femoral, Tibial New Axis , 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6660-01042	Femoral, Tibial New Axis , 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6660-01062	Femoral, Tibial New Axis , 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6660-01082	Femoral, Tibial New Axis , 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6660-01102	Femoral, Tibial New Axis , 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6663-01021	Tibia, Total Condyle Axis 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6663-01041	Tibia, Total Condyle Axis 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6663-01061	Tibia, Total Condyle Axis 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6663-01081	Tibia, Total Condyle Axis 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6663-01101	Tibia, Total Condyle Axis 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6663-01022	Tibia, Total Condyle Axis 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6663-01042	Tibia, Total Condyle Axis 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6663-01062	Tibia, Total Condyle Axis 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6663-01082	Tibia, Total Condyle Axis 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6663-01102	Tibia, Total Condyle Axis 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6665-01021	Femoral New Axis , 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6665-01041	Femoral New Axis , 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6665-01061	Femoral New Axis , 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6665-01081	Femoral New Axis , 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE

6665-01101	Femoral New Axis , 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6665-01022	Femoral New Axis , 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6665-01042	Femoral New Axis , 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6665-01062	Femoral New Axis , 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6665-01082	Femoral New Axis , 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6665-01102	Femoral New Axis , 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6666-01021	Tibial New Axis , 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6666-01041	Tibial New Axis , 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6666-01061	Tibial New Axis , 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6666-01081	Tibial New Axis , 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6666-01101	Tibial New Axis , 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6666-01022	Tibial New Axis , 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6666-01042	Tibial New Axis , 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6666-01062	Tibial New Axis , 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6666-01082	Tibial New Axis , 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6666-01102	Tibial New Axis , 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6681-01021	Total Femoral Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6681-01022	Total Femoral Axis, 2R, Condyle Width 58 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6681-01041	Total Femoral Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6681-01042	Total Femoral Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6681-01061	Total Femoral Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6681-01062	Total Femoral Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6681-01081	Total Femoral Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6681-01082	Total Femoral Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6681-01101	Total Femoral Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE

6681-01102	Total Femoral Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6690-01090	Custom Centralizer, 9#, Diameter 9 mm	9	UHMWPE
6690-01100	Custom Centralizer, 10#, Diameter 10 mm	10	UHMWPE
6690-01110	Custom Centralizer, 11#, Diameter 11 mm	11	UHMWPE
6690-01120	Custom Centralizer, 12#, Diameter 12 mm	12	UHMWPE
6690-01130	Custom Centralizer, 13#, Diameter 13 mm	13	UHMWPE
6690-01140	Custom Centralizer, 14#, Diameter 14 mm	14	UHMWPE
6690-01150	Custom Centralizer, 15#, Diameter 15 mm	15	UHMWPE



Knee Joint Prosthesis-New Axial (Custom Ordinary)

Code	Description	width	Material
6640-02021	Femoral Adjustable New Axis 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6640-02041	Femoral Adjustable New Axis 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6640-02061	Femoral Adjustable New Axis 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6640-02081	Femoral Adjustable New Axis 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6640-02101	Femoral Adjustable New Axis 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6640-02022	Femoral Adjustable New Axis 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6640-02042	Femoral Adjustable New Axis 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6640-02062	Femoral Adjustable New Axis 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6640-02082	Femoral Adjustable New Axis 8R, Condyle Width 71 mm, Tray width 77 mm	66, 72	Titanium, CoCrMo, UHMWPE
6640-02102	Femoral Adjustable New Axis 10R, Condyle Width 76 mm, Tray width 82 mm	66, 72	Titanium, CoCrMo, UHMWPE
6642-02021	Tibial Adjustable New Axis 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6642-02041	Tibial Adjustable New Axis 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6642-02061	Tibial Adjustable New Axis 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6642-02081	Tibial Adjustable New Axis 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE

6642-02101	Tibial Adjustable New Axis 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6642-02022	Tibial Adjustable New Axis 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6642-02042	Tibial Adjustable New Axis 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6642-02062	Tibial Adjustable New Axis 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6642-02082	Tibial Adjustable New Axis 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6642-02102	Tibial Adjustable New Axis 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6644-02021	Child Reserved Epiphyseal Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6644-02041	Child Reserved Epiphyseal Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6644-02061	Child Reserved Epiphyseal Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6644-02081	Child Reserved Epiphyseal Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6644-02101	Child Reserved Epiphyseal Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6644-02022	Child Reserved Epiphyseal Axis, 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6644-02042	Child Reserved Epiphyseal Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6644-02062	Child Reserved Epiphyseal Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6644-02082	Child Reserved Epiphyseal Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6644-02102	Child Reserved Epiphyseal Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6645-02021	Femoral Adjustable Half Joint, 2L, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6645-02041	Femoral Adjustable Half Joint, 4L, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6645-02061	Femoral Adjustable Half Joint, 6L, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6645-02081	Femoral Adjustable Half Joint, 8L, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6645-02101	Femoral Adjustable Half Joint, 10L, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6645-02022	Femoral Adjustable Half Joint, 2R, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6645-02042	Femoral Adjustable Half Joint, 4R, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6645-02062	Femoral Adjustable Half Joint, 6R, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6645-02082	Femoral Adjustable Half Joint, 8R, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE

6645-02102	Femoral Adjustable Half Joint, 10R, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6646-02021	Tibial Adjustable Half Joint, 2LTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6646-02041	Tibial Adjustable Half Joint, 4LTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6646-02061	Tibial Adjustable Half Joint, 6LTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6646-02081	Tibial Adjustable Half Joint, 8LTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6646-02101	Tibial Adjustable Half Joint, 10LTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6646-02022	Tibial Adjustable Half Joint, 2RTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6646-02042	Tibial Adjustable Half Joint, 4RTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6646-02062	Tibial Adjustable Half Joint, 6RTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6646-02082	Tibial Adjustable Half Joint, 8RTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6646-02102	Tibial Adjustable Half Joint, 10RTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6647-02021	Tibia, Total Femoral Adjustable New Axis, 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6647-02041	Tibia, Total Femoral Adjustable New Axis, 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6647-02061	Tibia, Total Femoral Adjustable New Axis, 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6647-02081	Tibia, Total Femoral Adjustable New Axis, 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6647-02101	Tibia, Total Femoral Adjustable New Axis, 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6647-02022	Tibia, Total Femoral Adjustable New Axis, 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6647-02042	Tibia, Total Femoral Adjustable New Axis, 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6647-02062	Tibia, Total Femoral Adjustable New Axis, 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6647-02082	Tibia, Total Femoral Adjustable New Axis, 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6647-02102	Tibia, Total Femoral Adjustable New Axis, 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6648-02021	Femoral Epiphysis Reserved Adjustable Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6648-02041	Femoral Epiphysis Reserved Adjustable Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6648-02061	Femoral Epiphysis Reserved Adjustable Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6648-02081	Femoral Epiphysis Reserved Adjustable Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE

6648-02101	Femoral Epiphysis Reserved Adjustable Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6648-02022	Femoral Epiphysis Reserved Adjustable Axis, 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6648-02042	Femoral Epiphysis Reserved Adjustable Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6648-02062	Femoral Epiphysis Reserved Adjustable Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6648-02082	Femoral Epiphysis Reserved Adjustable Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6648-02102	Femoral Epiphysis Reserved Adjustable Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6658-02021	Half Femoral Joint, 2L, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6658-02041	Half Femoral Joint, 4L, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6658-02061	Half Femoral Joint, 6L, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6658-02081	Half Femoral Joint, 8L, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6658-02101	Half Femoral Joint, 10L, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6658-02022	Half Femoral Joint, 2R, Condyle Width 56 mm	56	Titanium, CoCrMo, UHMWPE
6658-02042	Half Femoral Joint, 4R, Condyle Width 61 mm	61	Titanium, CoCrMo, UHMWPE
6658-02062	Half Femoral Joint, 6R, Condyle Width 66 mm	66	Titanium, CoCrMo, UHMWPE
6658-02082	Half Femoral Joint, 8R, Condyle Width 71 mm	71	Titanium, CoCrMo, UHMWPE
6658-02102	Half Femoral Joint, 10R, Condyle Width 76 mm	76	Titanium, CoCrMo, UHMWPE
6659-02021	Half Tibial Joint, 2LTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6659-02041	Half Tibial Joint, 4LTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6659-02061	Half Tibial Joint, 6LTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6659-02081	Half Tibial Joint, 8LTray width 72 mm	72	Titanium, CoCrMo, UHMWPE
6659-02101	Half Tibial Joint, 10LTray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6659-02022	Half Tibial Joint, 2RTray width 58 mm	58	Titanium, CoCrMo, UHMWPE
6659-02042	Half Tibial Joint, 4RTray width 62 mm	62	Titanium, CoCrMo, UHMWPE
6659-02062	Half Tibial Joint, 6RTray width 67 mm	67	Titanium, CoCrMo, UHMWPE
6659-02082	Half Tibial Joint, 8RTray width 72 mm	72	Titanium, CoCrMo, UHMWPE

6659-02102	Half Tibial Joint, 10R Tray width 77 mm	77	Titanium, CoCrMo, UHMWPE
6660-02021	Femoral, Tibial New Axis , 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6660-02041	Femoral, Tibial New Axis , 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6660-02061	Femoral, Tibial New Axis , 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6660-02081	Femoral, Tibial New Axis , 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6660-02101	Femoral, Tibial New Axis , 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6660-02022	Femoral, Tibial New Axis , 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6660-02042	Femoral, Tibial New Axis , 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6660-02062	Femoral, Tibial New Axis , 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6660-02082	Femoral, Tibial New Axis , 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6660-02102	Femoral, Tibial New Axis , 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6663-02021	Tibia, Total Condyle Axis 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6663-02041	Tibia, Total Condyle Axis 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6663-02061	Tibia, Total Condyle Axis 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6663-02081	Tibia, Total Condyle Axis 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6663-02101	Tibia, Total Condyle Axis 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6663-02022	Tibia, Total Condyle Axis 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6663-02042	Tibia, Total Condyle Axis 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6663-02062	Tibia, Total Condyle Axis 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6663-02081	Tibia, Total Condyle Axis 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6663-02101	Tibia, Total Condyle Axis 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6665-02021	Femoral New Axis , 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6665-02041	Femoral New Axis , 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6665-02061	Femoral New Axis , 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6665-02081	Femoral New Axis , 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE

6665-02101	Femoral New Axis , 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6665-02022	Femoral New Axis , 2R, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6665-02042	Femoral New Axis , 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6665-02062	Femoral New Axis , 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6665-02082	Femoral New Axis , 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6665-02102	Femoral New Axis , 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6666-02021	Tibial New Axis , 2L, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6666-02041	Tibial New Axis , 4L, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6666-02061	Tibial New Axis , 6L, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6666-02081	Tibial New Axis , 8L, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6666-02101	Tibial New Axis , 10L, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6666-02022	Tibial New Axis , 2R, Condyle Width 56 mm, Tray width 58 mm	56, 58	Titanium, CoCrMo, UHMWPE
6666-02042	Tibial New Axis , 4R, Condyle Width 61 mm, Tray width 62 mm	61, 62	Titanium, CoCrMo, UHMWPE
6666-02062	Tibial New Axis , 6R, Condyle Width 66 mm, Tray width 67 mm	66, 67	Titanium, CoCrMo, UHMWPE
6666-02082	Tibial New Axis , 8R, Condyle Width 71 mm, Tray width 72 mm	71, 72	Titanium, CoCrMo, UHMWPE
6666-02102	Tibial New Axis , 10R, Condyle Width 76 mm, Tray width 77 mm	76, 77	Titanium, CoCrMo, UHMWPE
6681-02021	Total Femoral Axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6681-02022	Total Femoral Axis, 2R, Condyle Width 58 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6681-02041	Total Femoral Axis, 4L, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6681-02042	Total Femoral Axis, 4R, Condyle Width 61 mm, Tray width 67 mm	61, 67	Titanium, CoCrMo, UHMWPE
6681-02061	Total Femoral Axis, 6L, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6681-02062	Total Femoral Axis, 6R, Condyle Width 66 mm, Tray width 72 mm	66, 72	Titanium, CoCrMo, UHMWPE
6681-02081	Total Femoral Axis, 8L, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6681-02082	Total Femoral Axis, 8R, Condyle Width 71 mm, Tray width 77 mm	71, 77	Titanium, CoCrMo, UHMWPE
6681-02101	Total Femoral Axis, 10L, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE

6681-02102	Total Femoral Axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
6690-02090	Custom Centralizer, 9#, Diameter 9 mm	9	UHMWPE
6690-02100	Custom Centralizer, 10#, Diameter 10 mm	10	UHMWPE
6690-02110	Custom Centralizer, 11#, Diameter 11 mm	11	UHMWPE
6690-02120	Custom Centralizer, 12#, Diameter 12 mm	12	UHMWPE
6690-02130	Custom Centralizer, 13#, Diameter 13 mm	13	UHMWPE
6690-02140	Custom Centralizer, 14#, Diameter 14 mm	14	UHMWPE
6690-02150	Custom Centralizer, 15#, Diameter 15 mm	15	UHMWPE



Knee Joint Prosthesis-Simple Knee Joint (Revision)

Code	Description	width	Material
6656-01021	Pure new axis, 2L, Condyle Width 56 mm, Tray width 62 mm	56, 62	Titanium, CoCrMo, UHMWPE
6656-01041	Pure new axis, 4L, Condyle Width 61 mm, Tray width 67 mm	56, 62	Titanium, CoCrMo, UHMWPE
6656-01061	Pure new axis, 6L, Condyle Width 66 mm, Tray width 72 mm	61, 67	Titanium, CoCrMo, UHMWPE
6656-01081	Pure new axis, 8L, Condyle Width 71 mm, Tray width 77 mm	61, 67	Titanium, CoCrMo, UHMWPE
6656-01101	Pure new axis, 10L, Condyle Width 76 mm, Tray width 82 mm	66, 72	Titanium, CoCrMo, UHMWPE
6656-01022	Pure new axis, 2R, Condyle Width 56 mm, Tray width 62 mm	66, 72	Titanium, CoCrMo, UHMWPE
6656-01042	Pure new axis, 4R, Condyle Width 61 mm, Tray width 67 mm	71, 77	Titanium, CoCrMo, UHMWPE
6656-01062	Pure new axis, 6R, Condyle Width 66 mm, Tray width 72 mm	71, 77	Titanium, CoCrMo, UHMWPE
6656-01082	Pure new axis, 8R, Condyle Width 71 mm, Tray width 77 mm	76, 82	Titanium, CoCrMo, UHMWPE
6656-01102	Pure new axis, 10R, Condyle Width 76 mm, Tray width 82 mm	76, 82	Titanium, CoCrMo, UHMWPE
New hip joint-upper hip			
3341-02090	New upper hip, Osteotomy diameter 24 mm, Osteotomy Length 50 mm	24, 50	Titanium
Knee prosthesis -Osteotomy segment			
6272-24101	Osteotomy segment 24/100, Diameter 24 mm, Length 100 mm	100	Titanium,

6272-24111	Osteotomy segment 24/110, Diameter 24 mm, Length 110 mm	110	Titanium,
6272-24121	Osteotomy segment 24/120, Diameter 24 mm, Length 120 mm	120	Titanium,
6272-24041	Osteotomy segment 24/40, Diameter 24 mm, Length 40 mm	40	Titanium,
6272-24051	Osteotomy segment 24/50, Diameter 24 mm, Length 50 mm	50	Titanium,
6272-24061	Osteotomy segment 24/60, Diameter 24 mm, Length 60 mm	60	Titanium,
6272-24071	Osteotomy segment 24/70, Diameter 24 mm, Length 70 mm	70	Titanium,
6272-24081	Osteotomy segment 24/80, Diameter 24 mm, Length 80 mm	80	Titanium,
6272-24091	Osteotomy segment 24/90, Diameter 24 mm, Length 90 mm	90	Titanium,
Knee joint prosthesis-femoral stem			
6271-05311	Universal Stem 11/105/30Stem Diameter 11 mm, Stem Length 105 mm, With osteotomy segment 30 mm	105	Titanium,
6271-05411	Universal Stem 11/105/40Stem Diameter 11 mm, Stem Length 105 mm, With osteotomy segment 40 mm	105	Titanium,
6271-05511	Universal Stem 11/105/50Stem Diameter 11 mm, Stem Length 105 mm, With osteotomy segment 50 mm	105	Titanium,
6271-05611	Universal Stem 11/105/60Stem Diameter 11 mm, Stem Length 105 mm, With osteotomy segment 60 mm	105	Titanium,
6271-15311	Universal Stem 11/115/30Stem Diameter 11 mm, Stem Length 115 mm, With osteotomy segment 30 mm	115	Titanium,
6271-15411	Universal Stem 11/115/40Stem Diameter 11 mm, Stem Length 115 mm, With osteotomy segment 40 mm	115	Titanium,
6271-15511	Universal Stem 11/115/50Stem Diameter 11 mm, Stem Length 115 mm, With osteotomy segment 50 mm	115	Titanium,
6271-15611	Universal Stem 11/115/60Stem Diameter 11 mm, Stem Length 115 mm, With osteotomy segment 60 mm	115	Titanium,

6271-25311	Universal Stem 11/125/30Stem Diameter 11 mm, Stem Length 125 mm, With osteotomy segment 30 mm	125	Titanium,
6271-25411	Universal Stem 11/125/40Stem Diameter 11 mm, Stem Length 125 mm, With osteotomy segment 40 mm	125	Titanium,
6271-25511	Universal Stem 11/125/50Stem Diameter 11 mm, Stem Length 125 mm, With osteotomy segment 50 mm	125	Titanium,
6271-25611	Universal Stem 11/125/60Stem Diameter 11 mm, Stem Length 125 mm, With osteotomy segment 60 mm	125	Titanium,
6271-15321	Universal Stem 12/115/30Stem Diameter 12 mm, Stem Length 115 mm, With osteotomy segment 30 mm	115	Titanium,
6271-15421	Universal Stem 12/115/40Stem Diameter 12 mm, Stem Length 115 mm, With osteotomy segment 40 mm	115	Titanium,
6271-15521	Universal Stem 12/115/50Stem Diameter 12 mm, Stem Length 115 mm, With osteotomy segment 50 mm	115	Titanium,
6271-15621	Universal Stem 12/115/60Stem Diameter 12 mm, Stem Length 115 mm, With osteotomy segment 60 mm	115	Titanium,
6271-25321	Universal Stem 12/125/30Stem Diameter 12 mm, Stem Length 125 mm, With osteotomy segment 30 mm	125	Titanium,
6271-25421	Universal Stem 12/125/40Stem Diameter 12 mm, Stem Length 125 mm, With osteotomy segment 40 mm	125	Titanium,
6271-25521	Universal Stem 12/125/50Stem Diameter 12 mm, Stem Length 125 mm, With osteotomy segment 50 mm	125	Titanium,
6271-25621	Universal Stem 12/125/60Stem Diameter 12 mm, Stem Length 125 mm, With osteotomy segment 60 mm	125	Titanium,
6271-35321	Universal Stem 12/135/30Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,

6271-35421	Universal Stem 12/135/40Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35521	Universal Stem 12/135/50Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35621	Universal Stem 12/135/60Stem Diameter 12 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-25331	Universal Stem 13/125/30Stem Diameter 13 mm, Stem Length 125 mm, With osteotomy segment 30 mm	125	Titanium,
6271-25431	Universal Stem 13/125/40Stem Diameter 13 mm, Stem Length 125 mm, With osteotomy segment 40 mm	125	Titanium,
6271-25531	Universal Stem 13/125/50Stem Diameter 13 mm, Stem Length 125 mm, With osteotomy segment 50 mm	125	Titanium,
6271-25631	Universal Stem 13/125/60Stem Diameter 13 mm, Stem Length 125 mm, With osteotomy segment 60 mm	125	Titanium,
6271-35331	Universal Stem 13/135/30Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 30 mm	135	Titanium,
6271-35431	Universal Stem 13/135/40Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 40 mm	135	Titanium,
6271-35531	Universal Stem 13/135/50Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 50 mm	135	Titanium,
6271-35631	Universal Stem 13/135/60Stem Diameter 13 mm, Stem Length 135 mm, With osteotomy segment 60 mm	135	Titanium,
6271-45331	Universal Stem 13/145/30Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,
6271-45431	Universal Stem 13/145/40Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 30 mm	145	Titanium,

6271-45531	Universal Stem 13/145/50Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 50 mm	145	Titanium,
6271-45631	Universal Stem 13/145/60Stem Diameter 13 mm, Stem Length 145 mm, With osteotomy segment 60 mm	145	Titanium,



XK Knee System



Femoral Condyle

Code	Description	width	Material
5310-04020	Cemented femoral condyle (CLGMC), 2#, width 17.5mm	17.5	CoCrMo
5310-04040	Cemented femoral condyle (CLGMC), 4#, width 19mm	19	CoCrMo
5310-04060	Cemented femoral condyle (CLGMC), 6#, width 20.5mm	20.5	CoCrMo
5310-04080	Cemented femoral condyle (CLGMC), 8#, width 22mm	22	CoCrMo



Tibial tray

Code	Description	width	Material
5311-04021	Cemented tibial tray (CLJMC), 2#L, width 24mm	24	CoCrMo
5311-04022	Cemented tibial tray (CLJMC), 2#R, width 24mm	24	CoCrMo
5311-04031	Cemented tibial tray (CLJMC), 3#L, width 26mm	26	CoCrMo
5311-04032	Cemented tibial tray (CLJMC), 3#R, width 26mm	26	CoCrMo
5311-04041	Cemented tibial tray (CLJMC), 4#L, width 26.5mm	26.5	CoCrMo
5311-04042	Cemented tibial tray (CLJMC), 4#R, width 26.5mm	26.5	CoCrMo
5311-04051	Cemented tibial tray (CLJMC), 5#L, width 28mm	28	CoCrMo
5311-04052	Cemented tibial tray (CLJMC), 5#R, width 28mm	28	CoCrMo
5311-04061	Cemented tibial tray (CLJMC), 6#L, width 30mm	30	CoCrMo
5311-04062	Cemented tibial tray (CLJMC), 6#R, width 30mm	30	CoCrMo
5311-04071	Cemented tibial tray (CLJMC), 7#L, width 31.5mm	31.5	CoCrMo
5311-04072	Cemented tibial tray (CLJMC), 7#R, width 31.5mm	31.5	CoCrMo
5311-04081	Cemented tibial tray (CLJMC), 8#L, width 33.5mm	33.5	CoCrMo
5311-04082	Cemented tibial tray (CLJMC), 8#R, width 33.5mm	33.5	CoCrMo



Tibial insert

Code	Description	width	Material
5312-03021	Tibial insert (CLJMP), 2#L, thickness 3mm, width 20.5mm	20.5	UHMWPE
5312-03022	Tibial insert (CLJMP), 2#R, thickness 3mm, width 20.5mm	20.5	UHMWPE
5312-04021	Tibial insert (CLJMP), 2#L, thickness 4mm, width 20.5mm	20.5	UHMWPE
5312-04022	Tibial insert (CLJMP), 2#R, thickness 4mm, width 20.5mm	20.5	UHMWPE
5312-05021	Tibial insert (CLJMP), 2#L, thickness 5mm, width 20.5mm	20.5	UHMWPE
5312-05022	Tibial insert (CLJMP), 2#R, thickness 5mm, width 20.5mm	20.5	UHMWPE
5312-06021	Tibial insert (CLJMP), 2#L, thickness 6mm, width 20.5mm	20.5	UHMWPE
5312-06022	Tibial insert (CLJMP), 2#R, thickness 6mm, width 20.5mm	20.5	UHMWPE
5312-07021	Tibial insert (CLJMP), 2#L, thickness 7mm, width 20.5mm	20.5	UHMWPE
5312-07022	Tibial insert (CLJMP), 2#R, thickness 7mm, width 20.5mm	20.5	UHMWPE
5312-08021	Tibial insert (CLJMP), 2#L, thickness 8mm, width 20.5mm	20.5	UHMWPE
5312-08022	Tibial insert (CLJMP), 2#R, thickness 8mm, width 20.5mm	20.5	UHMWPE
5312-09021	Tibial insert (CLJMP), 2#L, thickness 9mm, width 20.5mm	20.5	UHMWPE
5312-09022	Tibial insert (CLJMP), 2#R, thickness 9mm, width 20.5mm	20.5	UHMWPE
5312-03041	Tibial insert (CLJMP), 4#L, thickness 3mm, width 22mm	22	UHMWPE
5312-03042	Tibial insert (CLJMP), 4#R, thickness 3mm, width 22mm	22	UHMWPE

5312-04041	Tibial insert (CLJMP), 4#L, thickness 4mm, width 22mm	22	UHMWPE
5312-04042	Tibial insert (CLJMP), 4#R, thickness 4mm, width 22mm	22	UHMWPE
5312-05041	Tibial insert (CLJMP), 4#L, thickness 5mm, width 22mm	22	UHMWPE
5312-05042	Tibial insert (CLJMP), 4#R, thickness 5mm, width 22mm	22	UHMWPE
5312-06041	Tibial insert (CLJMP), 4#L, thickness 6mm, width 22mm	22	UHMWPE
5312-06042	Tibial insert (CLJMP), 4#R, thickness 6mm, width 22mm	22	UHMWPE
5312-07041	Tibial insert (CLJMP), 4#L, thickness 7mm, width 22mm	22	UHMWPE
5312-07042	Tibial insert (CLJMP), 4#R, thickness 7mm, width 22mm	22	UHMWPE
5312-08041	Tibial insert (CLJMP), 4#L, thickness 8mm, width 22mm	22	UHMWPE
5312-08042	Tibial insert (CLJMP), 4#R, thickness 8mm, width 22mm	22	UHMWPE
5312-09041	Tibial insert (CLJMP), 4#L, thickness 9mm, width 22mm	22	UHMWPE
5312-09042	Tibial insert (CLJMP), 4#R, thickness 9mm, width 22mm	22	UHMWPE
5312-03061	Tibial insert (CLJMP), 6#L, thickness 3mm, width 24mm	24	UHMWPE
5312-03062	Tibial insert (CLJMP), 6#R, thickness 3mm, width 24mm	24	UHMWPE
5312-04061	Tibial insert (CLJMP), 6#L, thickness 4mm, width 24mm	24	UHMWPE
5312-04062	Tibial insert (CLJMP), 6#R, thickness 4mm, width 24mm	24	UHMWPE
5312-05061	Tibial insert (CLJMP), 6#L, thickness 5mm, width 24mm	24	UHMWPE
5312-05062	Tibial insert (CLJMP), 6#R, thickness 5mm, width 24mm	24	UHMWPE
5312-06061	Tibial insert (CLJMP), 6#L, thickness 6mm, width 24mm	24	UHMWPE
5312-06062	Tibial insert (CLJMP), 6#R, thickness 6mm, width 24mm	24	UHMWPE
5312-07061	Tibial insert (CLJMP), 6#L, thickness 7mm, width 24mm	24	UHMWPE
5312-07062	Tibial insert (CLJMP), 6#R, thickness 7mm, width 24mm	24	UHMWPE
5312-08061	Tibial insert (CLJMP), 6#L, thickness 8mm, width 24mm	24	UHMWPE
5312-08062	Tibial insert (CLJMP), 6#R, thickness 8mm, width 24mm	24	UHMWPE
5312-09061	Tibial insert (CLJMP), 6#L, thickness 9mm, width 24mm	24	UHMWPE

5312-09062	Tibial insert (CLJMP), 6#R, thickness 9mm, width 24mm	24	UHMWPE
5312-03081	Tibial insert (CLJMP), 8#L, thickness 3mm, width 26mm	26	UHMWPE
5312-03082	Tibial insert (CLJMP), 8#R, thickness 3mm, width 26mm	26	UHMWPE
5312-04081	Tibial insert (CLJMP), 8#L, thickness 4mm, width 26mm	26	UHMWPE
5312-04082	Tibial insert (CLJMP), 8#R, thickness 4mm, width 26mm	26	UHMWPE
5312-05081	Tibial insert (CLJMP), 8#L, thickness 5mm, width 26mm	26	UHMWPE
5312-05082	Tibial insert (CLJMP), 8#R, thickness 5mm, width 26mm	26	UHMWPE
5312-06081	Tibial insert (CLJMP), 8#L, thickness 6mm, width 26mm	26	UHMWPE
5312-06082	Tibial insert (CLJMP), 8#R, thickness 6mm, width 26mm	26	UHMWPE
5312-07081	Tibial insert (CLJMP), 8#L, thickness 7mm, width 26mm	26	UHMWPE
5312-07082	Tibial insert (CLJMP), 8#R, thickness 7mm, width 26mm	26	UHMWPE
5312-08081	Tibial insert (CLJMP), 8#L, thickness 8mm, width 26mm	26	UHMWPE
5312-08082	Tibial insert (CLJMP), 8#R, thickness 8mm, width 26mm	26	UHMWPE
5312-09081	Tibial insert (CLJMP), 8#L, thickness 9mm, width 26mm	26	UHMWPE
5312-09082	Tibial insert (CLJMP), 8#R, thickness 9mm, width 26mm	26	UHMWPE